



## Accessories for 6.3026.210

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery 6.3026.210

Qty.	Order no.	Description	
1 pcs	6.1446.030	<b>Ball stopper / B-14/(15)</b> For the movable holding of buret tips in titration vessels	
1 pcs	6.1542.020	<b>Flat stopcock / PCTFE/PTFE</b>	
1 pcs	6.1543.060	<b>Tip / M6 thread</b>  Tip with M6 thread. Together with the 6.1805.160 capillary and the 6.1446.030 link stopper, it makes up the entire 6.1537.010 buret tip.	
1 pcs	6.1543.200	<b>Titration tip / M6 thread</b>  Titration tip with antidiffusion valve.	

---

**1 pcs**     **6.1563.030**     **Light protection for glass cylinder of 806 Exchange units**  
For 6.3026.XXX Exchange units



---

**1 pcs**     **6.1602.105**     **Bottle attachment / GL 45 / green**

For direct dosing from reagent bottles with GL 45 thread.



---

**1 pcs**     **6.1608.023**     **Amber glass bottle / 1000 mL / GL 45**  
For Exchange Units. Bottle for auxiliary solutions.



---

**1 pcs**     **6.1805.010**     **FEP tubing / M6 / 13 cm**  
With light and kink protection.



---

**1 pcs**     **6.1805.080**     **FEP tubing / M6 / 25 cm**  
With light and kink protection.



---

**1 pcs**    **6.1805.100**    **FEP tubing / M6 / 40 cm**  
With light and kink protection



---

**1 pcs**    **6.1819.020**    **Tube**  
For M6 threaded holes. For 6.1602.1XX bottle attachments



---

**1 pcs**    **6.2023.020**    **Clip for SGJ 14/15**  
Clip for SGJ 14/16



---

**1 pcs**    **6.2046.070**    **Card holder**  
Holder for marking labels on 6.3026.XXX Exchange Units.



## Optional accessories

Order no.	Description	
<b>6.1563.040</b> For 6.3026.XXX Exchange Units	<b>Thermostat jacket for 806 Exchange Units</b>	
<b>6.1576.210</b>	<b>Exchange unit base plate / 10 mL</b>	
<b>6.1608.030</b>	<b>Round glass bottle / 1000 mL / GL 45</b>	
<b>6.1608.040</b> For Exchange Units. Bottle for auxiliary solutions.	<b>PE bottle / 1000 mL / GL 45</b>	