

Accessories for 2.728.0034

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.728.0034

Qty.	Order no.	Description	
1 pcs	1.728.0010	728 Stirrer without stand For mounting on a stand rod 10 mm dia.; without base plate (or holder) and without stand rod; without power adapter, with fixed connecting cable.	
1 pcs	6.1903.030	Stirring bar / 25 mm Stirring bar with magnetic core, PTFE covering, length 25 mm.	
1 pcs	6.2013.010	Clamping ring For support rods with a diameter of 10 mm.	
1 pcs	6.2026.010	Support stand with support rod Complete, with stand rod, total length 42 cm. Used with 649, 728 Magnetic Swing-out Stirrer alone (without Dosimat)	
1 pcs	6.2117.020	Power adapter EU 230 V / 9 V DC With EU mains plug. Input: 220...240 V AC. For 622, 649, 722, 728 Stirrers, 727 Ti Stand or Spectrosense. Output: 9 V DC/0.2 A	