




Accessories for 2.746.0020

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.746.0020

Qty.	Order no.	Description	
1 pcs	1.746.0020	746 VA Trace Analyzer with integrated printer The 746 VA Trace Analyzer is the compact controller for the 747 VA Stand and the 695 VA Autosampler.	
1 pcs	6.2125.110	RS-232 cable 3 m, 25-pin – 9-pin Connecting Metrohm RS-232 interface to PC. Used with pH Meters, Titrinos, KF Coulometers, VA Trace Analyzer.	
1 pcs	6.2237.040	Thermal paper 40 m Sufficient for approx. 10'000 printed lines.	

Optional accessories

Order no.	Description	
2.700.0010	700 Dosino Drive for Dosing Unit, with fixed cable and 9-pin plug (DB9).	
2.772.0110	772 Pump Unit Peristaltic pump for direct connection to Sample Processors.	
2.823.0010	823 Membrane Pump Unit Membrane pump for direct connection to Sample Processors.	
6.2134.000	Extension cable / 2 m / D-sub 9 pin (1:1) Connecting 693 VA Processor or 746 VA Trace Analyzer to 685 Dosimat or 762 IC Interface to 750 Autosampler 9-pin extension cable (1:1). Specially shielded.	
6.5327.000	MVA-Hg, equipment for mercury determination Complete set of accessories for the determination of mercury as per Application Bulletin 96. Contains rotating gold working electrode, reference electrode, glassy-carbon auxiliary electrode and other accessories.	