

## Accessories for 2.760.0010

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery 2.760.0010

Qty.	Order no.	Description	
2 pcs	6.1236.020	<b>Sleeve with SGJ 14/12 mm</b> Sleeve with SGJ 14/12 mm and O-ring.	
1 pcs	6.1446.160	<b>Thermometer stopper, B-14/15, oblique insertion</b> To have the inserts of the 6.1458.010 titration head in a sloping position.	
1 pcs	6.1805.110	<b>FEP tubing / M6 / 80 cm</b> With light and kink protection	
1 pcs	6.1815.010	<b>Spiral band 0.5 m</b> For holding together different cables or tubing.	

---

**1 pcs**     **6.2041.310**     **Sample rack 12 x 250 mL**  
Sample rack for 12 x 250 mL, one row



---

**3 pcs**     **6.2709.070**     **Guiding sleeve**  
For fixation of buret and dosing tips in titration heads and electrode holders



---

**1 pcs**     **6.2752.010**     **Plug cover**  
Used with the 717, 730, 760 Sample Changers, IC Sample Processors, 774 Oven Sample Processor



---

**1 pcs**     **6.5610.010**     **KFT equipment for 810/814 Sample Processor**

Karl Fischer equipment for 810 or 814 Sample Processor.



## Optional accessories

Order no.	Description	
<b>6.1453.220</b>	<b>Sample beaker, plastic (PP), 200 mL, 1 piece</b>	
<p>Sample beaker, plastic (polypropylene, PP), 200 mL, suitable for sample racks: 6.2041.310, 6.2041.400, 6.2041.430, 6.2041.510, 6.2041.820, 6.2068.020, and 6.2068.030, 1 piece</p>		
<b>6.1453.250</b>	<b>Sample beaker, plastic (PP), 250 mL, 1 piece</b>	
<p>Sample beaker, plastic (polypropylene, PP), 250 mL, suitable for sample racks: 6.2041.310, 6.2041.400, 6.2041.430, 6.2041.510, 6.2041.820, 6.2068.020, and 6.2068.030, 1 piece</p>		
<b>6.1458.010</b>	<b>Titration head, 6 x SGJ 14 and 3 x SGJ 9 openings</b>	
<p>With 6xSGJ 14 and 3x SGJ 9 apertures.</p>		
<b>6.1458.020</b>	<b>Titration head /4 x M10 connectors</b>	
<p>With 4 M10 screw nipples, including stoppers.</p>		
<b>6.1458.030</b>	<b>KF titration head /4 x M10 connectors</b>	
<p>With 4 M10 screw nipples, including stoppers.</p>		