

## Accessories for 2.809.3010

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery 2.809.3010

Qty.	Order no.	Description
1 pcs	1.802.0010	<b>802 Rod Stirrer for 727 Ti Stand</b>

Rod stirrer with separate 6.2137.000 Speed controller and 6.1909.010 Propeller stirrer.



1 pcs	1.804.0010	<b>804 Ti Stand without stand rod</b>
-------	------------	---------------------------------------

Titration stand and controller for 802 Stirrer. The 804 Ti Stand together with the optional 802 Stirrer provides an alternative to the magnetic stirrer. Without stand.



1 pcs	1.809.0010	<b>809 Titrande</b>
-------	------------	---------------------

High-end titrator for up to four 800 Dosino dispensing systems. Dynamic (DET) and monotonic (MET) titration, endpoint titration (SET), measurements with ion-selective electrodes (MEAS CONC), dosing functions with monitoring (DOS), liquid handling. With four MSB connections, one galvanically separate measuring interface, USB connection.



---

**1 pcs    6.0229.010    Solvotrode easyClean (fixed cable 1.2 m)**

Combined pH electrode with flexible easyClean diaphragm and fixed cable (1.2 m) for all non-aqueous acid/base titrations. The glass membrane is optimized for poorly conducting solutions, and thanks to the easy-to-clean easyClean diaphragm, the electrode is well suited for heavily contaminated samples, such as used oil. This electrode can be used with non-aqueous reference electrolytes (lithium chloride or tetraethylammonium bromide). Storage in the respective reference electrolyte.



---

**1 pcs    6.1546.030    Extractor rod**

For PTFE piston of dosing unit



---

**2 pcs    6.1608.023    Amber glass bottle / 1000 mL / GL 45**

For Exchange Units. Bottle for auxiliary solutions.



---

**1 pcs    6.1909.010    Stirring propeller / 96 mm**

Propeller stirrer, length from lower edge of SGJ: 96 mm. For use with beakers with 722, 802 Rod Stirrer



---

**1 pcs    6.2001.060    Support stand without support rod**

To attach an 801, 804, 803 Stirrer to a Titrando, Titrino plus, Dosimat plus or 856 and 867.



---

**1 pcs**      **6.2013.010**      **Clamping ring**

For support rods with a diameter of 10 mm.



---

**1 pcs**      **6.2016.070**      **Support rod / 400 mm**



---

**1 pcs**      **6.2021.020**      **Electrode holder**

Electrode holder for 4 electrodes and 2 buret tips



---

**2 pcs**      **6.2043.005**      **Holding clip for bottles**

Holding spring for reagent bottles in exchange units.



---

**1 pcs**      **6.2103.130**      **Adapter red 2 mm plug / 4 mm socket**

For connecting plug B (4 mm) to 2 mm socket.



---

**1 pcs**      **6.2103.140**      **Adapter black 2 mm plug / B socket 4 mm**  
For connecting plug B (4 mm) to 2 mm socket.



---

**1 pcs**      **6.2151.000**      **Cable USB A – mini-DIN 8-pin**  
Controller cable



---

**1 pcs**      **6.2320.000**      **TEABr 0.4 mol/L in ethylene glycol (250 mL)**  
Electrolyte solution TEABr (tetraethylammonium bromide in ethylene glycol),  $c(\text{TEABr}) = 0.4 \text{ mol/L}$



---

**1 pcs**      **6.2621.070**      **5 mm hex key for IC Sample Processors**



---

**1 pcs**      **6.2621.130**      **Hexagon key 2 mm**  
2 mm.



---

**2 pcs****6.3032.220****Dosing Unit 20 mL**

Dosing unit with integrated data chip with 20 mL glass cylinder and light protection, mountable on reagent bottle with ISO/DIN GL45 thread. FEP tubing connection, anti-diffusion buret tip.

