

## Accessories for 2.827.0015

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery 2.827.0015

Qty.	Order no.	Description	
1 pcs	6.2013.010	<b>Clamping ring</b>  For support rods with a diameter of 10 mm.	
1 pcs	6.2016.050	<b>Support rod / 300 mm</b>	
1 pcs	6.2021.020	<b>Electrode holder</b> Electrode holder for 4 electrodes and 2 buret tips	
1 pcs	6.2161.020	<b>Plug-in power supply unit US 120V/6V DC</b> Power supply for 827 pH lab	
1 pcs	6.2621.070	<b>5 mm hex key for IC Sample Processors</b>	