










## Accessories for 2.833.0050

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery 2.833.0050

Qty.	Order no.	Description	
1 pcs	6.1803.040	<b>PTFE capillary 0.5 mm i.d. / 1 m</b> Capillary for sample handling in IC.	
1 pcs	6.1803.050	<b>PTFE capillary 0.5 mm i.d. / 20 cm</b> Capillary for inline dialysis and filtration, for Dialysis Unit, IC Filtration Sample Processor	
2 pcs	6.1803.070	<b>PTFE capillary 0.97 mm i.d. / 40 cm</b> Capillary for filtration, for IC Filtration Sample Processor, IC Liquid Handling Filtration Unit	
2 pcs	6.1803.080	<b>PTFE capillary 0.97 mm i.d. / 1 m</b> Capillary for filtration, for IC Filtration Sample Processor, IC Liquid Handling Filtration Unit	

<b>1 pcs</b>	<b>6.1803.100</b>	<b>PTFE capillary, 0.5 mm ID / 40 cm</b> Capillary for filtration and dialysis	
<b>1 pcs</b>	<b>6.2128.180</b>	<b>Cable 762/830 - 752/753/754/833</b> Connecting IC Interface - IC Pump Unit, IC Suppressor Module, IC Dialysis Unit or IC Liquid Handling Unit	
<b>1 pcs</b>	<b>6.2143.200</b>	<b>Cable 732/819 - 752/754/833</b> Connecting IC Detector - Pump Unit or IC Detector - Dialysis Unit	
<b>1 pcs</b>	<b>6.2621.070</b>	<b>5 mm hex key for IC Sample Processors</b>	
<b>1 pcs</b>	<b>6.2714.020</b>	<b>Filtration membrane 1</b> Nominal pore size : 0.2 µm. Set of 50 pieces. For inline filtration	
<b>1 pcs</b>	<b>6.2729.110</b>	<b>Ultrafiltration cell</b> With 3 x 6.2744.060 Blind stoppers. For inline filtration	

---

**1 pcs**     **6.2744.000**     **Pressure screw PVDF**  
With UNF 10/32 connection. For IC instruments, inline dialysis (connection of PTFE capillaries)



---

**1 pcs**     **6.2744.010**     **Pressure screw 5x**  
With UNF 10/32 connection. For the connection of PEEK capillaries



---

**1 pcs**     **6.2744.030**     **Coupling olive/UNF 10/32, 4x**  
Connection between pressure screw and pump tubing. 4 pieces. For IC instruments with peristaltic pump





---

**1 pcs**     **6.2831.010**     **Plastic tweezers**  
For dialysis and filtration membranes



## Optional accessories

Order no.	Description	
<b>6.1803.030</b>	<b>PTFE capillary 0.5 mm i.d. / 3 m</b> Capillary for inline dialysis, for Dialysis Unit, IC Dialysis Sample Processor, IC Liquid Handling Dialysis Unit	
<b>6.1826.060</b>	<b>Pump tubing (orange/yellow)</b> For Compact IC , Personal IC , Basic IC	
<b>6.2755.000</b>	<b>Tubing cartridge</b> For the peristaltic pumps of IC instruments. Contact pressure adjustable.	