




Accessories for 2.843.0050

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.843.0050

| Qty. | Order no. | Description | |
|-------|------------|--|---|
| 1 pcs | 1.843.0020 | 843 Pump Station (membrane) The 843 Pump Station (membrane) has two built-in membrane pumps. These can be controlled directly via the interface using remote signals or manually by pressing a button. |  |
| 1 pcs | 6.1543.170 | Aspiration tip / M8 thread Aspiration tip for the Sample Processors. |  |
| 2 pcs | 6.1621.000 | Container 10 L As rinsing or waste container in automated systems. |  |

2 pcs **6.1621.100** **Cap for 10-L container**
Cap for 10-L container. For easy and safe transportation of filled containers.



4 pcs **6.1801.120** **PVC tubing / 4 mm / 6 mm / 2 m**



3 pcs **6.1805.060** **FEP tubing / M6 / 60 cm**
With light and kink protection



3 pcs **6.1805.420** **FEP tubing / M6 / 48 cm**
With light and kink protection



1 pcs **6.1805.510** **PTFE tubing / M8 / 60 cm**

With protection against kinks.



1 pcs 6.1808.170 Distributor for Sample Processors for rinsing and aspirating

Distributor for installation on Sample Processors for rinsing and aspirating when using externally connected pumps.



2 pcs 6.1812.000 PTFE tubing 4 mm / 6 mm / 4 m



2 pcs 6.1828.000 PVDF connection nipple

For 6.1621.000 canister.



1 pcs 6.2141.300 Remote cable 843 Pump Station to Sample Processor

Remote cable for the direct connection of the 843 Pump Station to the pump connectors of the Sample Processors




3 pcs 6.2740.020 Spray nozzle

For the fine-spraying of the rinsing solution.



Optional accessories

| Order no. | Description | |
|---|-----------------------|---|
| 6.2740.030 | Rinsing nozzle |  |
| <p>The rinsing solution is applied as a fine jet for the optimal removal of layers on the electrodes and titration accessories. For 6.5613.000 ECO titration equipment.</p> | | |
