

## Accessories for 2.861.0020

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery 2.861.0020

Qty.	Order no.	Description	
1 pcs	1.861.0020	<b>861 Advanced Compact IC with MSM II</b>  The Advanced Compact IC is suited for applications with or without suppression such as anion determinations as well as for methods combining conductivity with other types of detection.	
2 pcs	6.1602.150	<b>Bottle attachment / GL 45 - 3 x UNF 10/32</b>  For the connection of 1/16 in. capillaries. Used in conjunction with MSM auxiliary solutions and in Metrohm Inline Sample Preparation.	
1 pcs	6.1602.160	<b>Bottle attachment GL 45 for eluents and reagents</b>  For eluent, reagent, and rinsing bottles, with connectors for the adsorber tube and the aspiration tubing.	
2 pcs	6.1608.023	<b>Amber glass bottle / 1000 mL / GL 45</b>  For Exchange Units. Bottle for auxiliary solutions.	

---

**1 pcs**    **6.1608.070**    **Eluent bottle / 2 L / GL 45**



---

**1 pcs**    **6.1609.000**    **Adsorbing tube, large and bent**  
For filling with adsorber material.



---

**1 pcs**    **6.1803.020**    **PTFE capillary 0.97 mm i.d. / 5 m**  
For all IC instruments



---

**2 pcs**    **6.1816.020**    **Silicone tubing 6 mm i.d., 1 m**  
For waste tubing.



---

**1 pcs**    **6.1834.010**    **PTFE aspiration tube to IC Pump, 2.5 m**  
For Compact IC, Personal IC, Basic IC, 6.5324.000 bottle holder.



---

**1 pcs**     **6.2023.020**     **Clip for SGJ 14/15**  
Clip for SGJ 14/16



---

**1 pcs**     **6.2027.030**     **Column holder 8.5 mm**  
For mounting columns in column chamber of IC instruments.



---

**1 pcs**     **6.2027.040**     **Column holder 11.3 mm**  
For mounting columns in column chamber of IC instruments.



---

**1 pcs**     **6.2134.100**     **RS-232 cable, 5 m, DB9 – DB9**  
Connecting Metrohm RS-232 interface to Titrios, 756 KF Coulometer, 762 or 771 IC Interfaces, 790 Personal IC, 743 Rancimat, 763 PVC Thermomat, 873 Biodiesel Rancimat.



---

**1 pcs**     **6.2617.010**     **Tool for piston seal**  
For removing and assembling the piston seal for all standard pump heads



---

**1 pcs**     **6.2620.150**     **Pulsation absorber (metal-free)**  
For all ion chromatographs



---

**1 pcs**     **6.2621.000**     **Adjustable wrench**  
Max. opening: 20 mm. For IC instruments



---

**1 pcs**     **6.2621.030**     **Hexagon key 4 mm**



---

**3 pcs**     **6.2744.010**     **Pressure screw 5x**  
With UNF 10/32 connection. For the connection of PEEK capillaries



---

**1 pcs**     **6.2744.020**     **Luer/UNF coupling**  
For IC instruments



---

**1 pcs    6.2744.034    Coupling olive/UNF 10/32, 2x**

Connection between pressure screw and pump tubing. 2 pieces. For IC instruments with peristaltic pump.



---

**2 pcs    6.2744.180    Pump tubing connector with security lock and filter**

For the connection of a pump tubing and a capillary with built-in filter



---

**1 pcs    6.2816.020    Syringe 10 mL with Luer connection**

For various applications in IC and VA



---

**1 pcs    6.2821.090    Aspiration filter**

Pore size 20 µm, set of 5 pieces. For 6.1834.000 Aspiration tubing and 6.1821.040 and 6.1821.050 Filter tubes.



---

**1 pcs    6.2821.130    Spare filter for inline filter**

Spare filter for inline filter.

