





Accessories for 2.886.0110

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.886.0110

Qty.	Order no.	Description	
1 pcs	1.886.0010	Professional Thermostat / Reactor **	
1 pcs	6.1015.100	Metrosep BP 1 Guard/2.0 The Metrosep BP 1 Guard/2.0 is used to generate a sufficiently high working pressure in the flow path of post-column reagents. They are used in combination with a high-pressure pump for conveying the post-column reagent in the Professional Reactor.	
1 pcs	6.1803.150	PTFE Capillary i.d. 0.25 mm / 1 m	
2 pcs	6.1831.100	PEEK capillary 0.25 mm i.d., 1 m	

1 pcs **6.2151.020** **Cable USB A - USB B 1.8 m**
USB connecting cable



1 pcs **6.2621.100** **3 mm hex key for IC Sample Processors**



1 pcs **6.2621.140** **Hexagon key 2.5 mm**



1 pcs **6.2744.010** **Pressure screw 5x**
With UNF 10/32 connection. For the connection of PEEK capillaries



1 pcs **6.2744.014** **Pressure screw 2x**
With UNF 10/32 connection. For the connection of PEEK capillaries



1 pcs **6.2845.100** **Reactor plate complete to Professional Reactor**
Intelligent reactor plate with reactor for PCR applications.



Optional accessories

Order no.	Description	
6.1602.200	Bottle cap 2 x M6, 1 x M8 to 6.1608.120	
<p>Bottle cap for eluent bottles made of PP for eluents, post-column derivatization solutions, etc. The PP eluent bottle can withstand gas pressure when this bottle cap is used.</p>		
6.1608.120	PP eluent bottle, 2 L	
<p>Eluent bottle for applications where glass should be avoided.</p>		
6.1836.020	Preheating capillary with core, 1.44 m	
<p>Preheating capillary to Professional Reactor / Heater</p>		
6.2061.100	Bottle holder for Professional IC instruments	
<p>Instrument cover and bottle holder for Professional IC instruments.</p>		
6.2061.110	Tray with sensor for Professional IC Instruments	
6.2061.120	System Connector	
<p>Connects an IC module with the footprint of a Professional IC with a Compact IC</p>		

6.2744.330 **Y connector 3 x UNF 10/32**



6.2845.110 **HT reactor plate, complete, for Professional Reactor**

Reactor plate with reactor for PCR applications at high temperatures (> 120 °C). With PEEK reactor capillaries.



6.2845.200 **Reactor complete to 6.2845.100**

Spare reactor for the reactor plate



6.2845.600 **Column holder complete to Professional Thermostat**

Intelligent column holder for the Professional Thermostat



6.5333.000 **IQ/OQ Kit for IC**

The IQ/OQ Kit contains all parts and standard solutions required for IQ/OQ in ion chromatography.



6.6059.021 **MagIC Net 2.x and 3.x Multi: 1 additional license**

1 additional license



6.6059.022 **MagIC Net 2.x and 3.x Multi: 5 additional licenses**

5 additional licenses



6.6059.023 **MagIC Net 2.x and 3.x Multi: 10 additional licenses**

10 additional licenses

