

## Accessories for 2.928.0230

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

## Optional accessories

Order no.	Description	
6.7440.170	<b>NIRS SingleFiber optical fiber, 2 m</b>  One-channel optical fiber in 2 m length for NIRS Process Analyzers	
6.7440.480	<b>NIRS Process, I/O Communication box</b>  I/O Communication box for process instruments with 24 analog and 24 digital channels	
6.7450.070	<b>NIRS 99% reflection and wavelength standard process, 1 inch diameter</b>  NIRS 99% reflection and wavelength standard for process instrument with 1 inch probe diameter	
6.7450.080	<b>NIRS 99% reflection and wavelength standard for process instrument, 0.5 inch probe diameter</b>  NIRS 99% reflection and wavelength standard for process instrument, 0.5 inch probe diameter	
6.7450.090	<b>NIRS 99% reflection and wavelength standard for process instrument, 1 inch diameter, 45° bevel</b>  NIRS 99% reflection and wavelength standard for process instrument with 1 inch diameter and 45° angle	