




## Accessories for 6.01600.010

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

## Optional accessories

Order no.	Description	
<b>6.01602.000</b>	<b>Aspiration coupling to 6.01600.010</b>	
	Aspiration coupling as spare part to a 6.01600.010 OMNIS Liquid Adapter	
<b>6.02106.000</b>	<b>Liquid Adapter extension cable, 0.55 m</b>	
	Extension cable (0.55 m in length) for connecting an OMNIS Liquid Adapter to an OMNIS Titrator, Titration Module or Dosing Module.	
<b>6.02106.010</b>	<b>Liquid Adapter extension cable, 2 m</b>	
	Extension cable (2 m in length) for connecting an OMNIS Liquid Adapter to an OMNIS Titrator, Titration Module or Dosing Module.	