








Accessories for 6.0278.300

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.0278.300

Qty.	Order no.	Description	
1 pcs	6.1236.050	Sleeve with SGJ 14/12 mm Sleeve with SGJ 14/12 mm, polyethylene.	
1 pcs	6.2008.040	Storage vessel Together with 6.2008.050 storage vessel holder provides a support for the electrode on 6.3032.XXX dosing units.	

Optional accessories

Order no.	Description	
2.854.0010	854 iConnect	
	854 iConnect - electrode cable and measuring amplifier for intelligent electrodes «iTrodes».	
6.2308.020	Electrolyte 3 mol/L KCl (250 mL)	
	Electrolyte solution c(KCl) = 3 mol/L, (for Ag/AgCl reference systems)	
6.2308.040	Idrolyte (250 mL)	
	Idrolyte (for 6.0224.100 Electrode and Unitrode at T >80°C)	
6.2323.000	Storage solution	
	Storage solution for all combined pH glass electrodes with reference electrolyte c(KCl) = 3 mol/L	
6.2325.000	pHit kit	
	Maintenance kit for pH electrodes	
	The kit contains:	
	<ul style="list-style-type: none"> • 50 mL cleaning solution • 50 mL 3M KCl solution • 50 mL storage solution • 2 Storage vessels • Instructions for use 	

6.2325.100**Cleaning solution 3 x 50 mL**

Reliable measuring results over long periods of time can be guaranteed only if the pH glass membrane and the diaphragm receive preventive and regular care. Cleaning by etching with toxic chemicals or applying mechanical treatment to the diaphragm is not only complicated and expensive, it also accelerates the ageing of the pH glass electrode.

The cleaning solution was developed for easy and gentle cleaning of pH glass electrodes. Regular use can considerably prolong their service life.

This cleaning solution is also part of the pHit kit.

