



Accessories for 6.1005.320

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Optional accessories

Order no.	Description	
6.1011.030	Metrosep RP 2 Guard/3.5 General guard column. For Metrosep Anion Dual 2 and Dual 3 in particular.	
6.1011.040	Metrosep RP 3 Guard HC/4.0 The Metrosep RP 3 Guard HC/4.0 is a guard column for universal use. It reliably protects the analytical separation column against contamination, securely removing lipophilic organic contamination, the smallest particles, traces of iron oxide and bacteria. The guard column is based on a polymer material and thanks to its larger pack volume, has a significantly higher capacity than the Metrosep RP 2Guard/3.5. It can be used in the entire pH range.	
6.1011.130	Spare filter for RP 2 Guard/3.5, 10 pieces	