

Accessories for 6.1050.420

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.



Scope of delivery 6.1050.420

Qty.	Order no.	Description
2 pcs	6.2744.060	Threaded stopper

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



Optional accessories

Order no.	Description
6.1050.500	Metrosep C 4 Guard/4.0
<p>The Metrosep C 4 Guard/4.0 contains the C-4 column material and is used to protect all Metrosep cation columns that have a substrate based on silica gel. Particles and contaminations are reliably retained, allowing the lifetime of the analytical separation column to be prolonged considerably. The economical price is an additional plus. The Metrosep C 4 Guard/4.0 also functions according to the "On Column Guard System" and is screwed directly onto the separation column with nearly no dead volume.</p>	
6.1050.510	Metrosep C 4 S-Guard/4.0
<p>The Metrosep C 4 S-Guard/4.0 contains the C-4 column material and is used to protect all Metrosep cation columns that have a substrate based on silica gel. Particles and contaminations are reliably retained, allowing the lifetime of the analytical separation column to be prolonged considerably. The economical price is an additional plus.</p>	
6.1050.530	Metrosep C 4 S-Guard - 50/4.0
<p>The Metrosep C 4 S-Guard - 50/4.0 contains the C 4 column material and is used to protect all Metrosep cation columns that have a substrate based on silica gel. Particles and contaminations are reliably retained, allowing the service life of the analytical separation column to be prolonged considerably. In comparison with the Standard C 4 guard columns, the Metrosep C 4 S-Guard - 50/4.0 exhibits greater capacity and therefore an even longer service life.</p>	