







Accessories for 6.1256.030

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.1256.030

Qty.	Order no.	Description	
1 pcs	6.1829.030	PTFE tubing / M6 / 21 cm, tapered	
1 pcs	6.1911.000	Flexible drive shaft for stirrer 10 pcs, for Professional VA instruments, including wrench.	
1 pcs	6.2069.000	Park station for the measuring head of Professional VA/CVS instruments	
1 pcs	6.2709.110	Stopper for the electrode opening on Professional VA/CVS instruments Set with 1 stopper for electrode opening (e.g. for pH electrode or Pt1000 temperature sensor), 6 stoppers for openings (auxiliary solutions and/or automatic sample addition).	

Optional accessories

Order no.	Description	
6.1454.130	Form seal made of FPM for 61414xxx	
	for KF and VA instruments	
6.2709.130	Adapter B-14/15 - SPE measuring head insert (VA)	
	Adapter for use of the electrode shaft (scTRACE Gold) 6.1241.080 in the 884 Professional VA with SPE measuring head 6.1256.030 or in the 797 VA Computrace or 663 VA Stand with measuring head insert 6.1412.020.	