

Accessories for 6.1462.160

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Optional accessories

Order no.	Description	
6.1458.010	Titration head, 6 x SGJ 14 and 3 x SGJ 9 openings With 6xSGJ 14 and 3x SGJ 9 apertures.	
6.1458.020	Titration head /4 x M10 connectors With 4 M10 screw nipples, including stoppers.	
6.1458.030	KF titration head /4 x M10 connectors With 4 M10 screw nipples, including stoppers.	
6.1458.040	Titration Head, 3 x SGJ 14 Titration head for 120 mL beakers. With 3 openings SGJ 14 and 4 openings 6.4 mm, 2 dosing tips included. Suitable for rinsing and aspiration equipment.	
6.1458.070	Titration head / 2 x SGJ 14 Titration head for 75 mL beakers. With 2 openings SGJ 14, 2 openings 6.4 mm, 1 opening M10 and 3 pre-mounted dosing tips. Suitable for rinsing and aspiration equipment.	