

Accessories for 6.5000.010

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.5000.010

Qty.	Order no.	Description	
1 pcs	6.1602.160	Bottle attachment GL 45 for eluents and reagents For eluent, reagent, and rinsing bottles, with connectors for the adsorber tube and the aspiration tubing.	
1 pcs	6.1608.070	Eluent bottle / 2 L / GL 45	
1 pcs	6.1609.000	Adsorbing tube, large and bent For filling with adsorber material.	
1 pcs	6.1803.040	PTFE capillary 0.5 mm i.d. / 1 m Capillary for sample handling in IC.	

1 pcs **6.2023.020** **Clip for SGJ 14/15**
Clip for SGJ 14/16



1 pcs **6.2744.014** **Pressure screw 2x**
With UNF 10/32 connection. For the connection of PEEK capillaries



1 pcs **6.2744.040** **2 x UNF 10/32 coupling**
For connecting 1/16 in. capillaries. For IC instruments



1 pcs **6.2744.070** **Pressure screw short**
Short version. With UNF 10/32 connection. 5 pieces. For the connection of PEEK capillaries



1 pcs **6.2744.090** **Pressure screw long**
Long version. With UNF 10/32 connection. 2 pieces. For the connection of PEEK capillaries (MCS and Sample degasser)



1 pcs **6.2744.210** **Tubing adapter for aspiration filter (ProfIC)**
For Professional IC instruments



1 pcs **6.2816.020** **Syringe 10 mL with Luer connection**
For various applications in IC and VA



1 pcs **6.2816.040** **Purge tubing**
With PTFE tubing and Luer connection. For syringes. For aspirating eluents.



1 pcs **6.2821.090** **Aspiration filter**
Pore size 20 µm, set of 5 pieces. For 6.1834.000 Aspiration tubing and 6.1821.040 and 6.1821.050 Filter tubes.



1 pcs **6.2821.130** **Spare filter for inline filter**
Spare filter for inline filter.

