


Accessories for 6.5000.020

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.5000.020

| Qty. | Order no. | Description | |
|-------|------------|--|---|
| 1 pcs | 6.1602.150 | Bottle attachment / GL 45 - 3 x UNF 10/32 |  |
| | | For the connection of 1/16 in. capillaries. Used in conjunction with MSM auxiliary solutions and in Metrohm Inline Sample Preparation. | |
| 1 pcs | 6.1608.020 | Glass bottle / 1000 mL / GL 45 |  |
| | | Bottle for auxiliary solutions. | |
| 1 pcs | 6.1803.030 | PTFE capillary 0.5 mm i.d. / 3 m |  |
| | | Capillary for inline dialysis, for Dialysis Unit, IC Dialysis Sample Processor, IC Liquid Handling Dialysis Unit | |
| 1 pcs | 6.2057.080 | Adsorption cartridge holder (ProfIC) |  |
| | | Holder to install adsorption cartridges in Professional IC Instruments | |

2 pcs **6.2744.040** **2 x UNF 10/32 coupling**
For connecting 1/16 in. capillaries. For IC instruments



1 pcs **6.2744.070** **Pressure screw short**
Short version. With UNF 10/32 connection. 5 pieces. For the connection of PEEK capillaries



1 pcs **6.2744.090** **Pressure screw long**
Long version. With UNF 10/32 connection. 2 pieces. For the connection of PEEK capillaries (MCS and Sample degasser)



1 pcs **6.2835.010** **Connecting piece to MSM and SPM module**

Universal connecting piece for MSM in Professional IC Vario and Compact Flex instruments and for MSM-HC and SPM module. With inlets and outlets.



1 pcs **6.2837.100** **CO2 adsorption cartridge CW**

Adsorption cartridge for removing the CO₂ from the air aspirated by the MCS. The cartridge is used without upstream H₂O adsorption cartridge.

