

Accessories for 6.5331.130

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.5331.130

Qty.	Order no.	Description	
1 pcs	6.1254.130	Platinum Variocell for Bioscan Platinum measuring cell with replaceable working electrode.	
1 pcs	6.1254.410	Spacer for Variocell, 50 µm 50 µm Spacer for 6.1254.1X0 Variocell	
1 pcs	6.1803.030	PTFE capillary 0.5 mm i.d. / 3 m Capillary for inline dialysis, for Dialysis Unit, IC Dialysis Sample Processor, IC Liquid Handling Dialysis Unit	
1 pcs	6.2744.014	Pressure screw 2x With UNF 10/32 connection. For the connection of PEEK capillaries	