




Accessories for 6.5339.070

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

Scope of delivery 6.5339.070

Qty.	Order no.	Description	
1 pcs	6.1204.500	Stirrer for Professional VA instruments	
		Stirrer for Professional VA instruments, complete, metal-free.	
1 pcs	6.1241.090	Electrode shaft for screen-printed electrodes (SPE)	
		Holder for <i>screen-printed electrodes</i> (SPE) with plug contacts for connecting the electrode cables for working, reference, and auxiliary electrodes. When used with screen-printed electrodes, the electrode shaft provides a complete sensor for voltammetry. Whereas the screen-printed electrodes must be replaced at regular intervals, the electrode shaft can continue to be used.	
1 pcs	6.1244.060	Silicone seal for SPE electrode shaft (10 pieces)	
		Seal made of silicone for SPE electrode shaft 6.1241.090 and measuring head (SPE) 6.01256.020, 10 pieces.	
1 pcs	6.1415.150	Titration vessel / 5-70 mL	