





Accessories for 2.889.0010

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.889.0010

Qty.	Order no.	Description	
1 pcs	1.889.0010	889 IC Sample Center	 <p>The 889 IC Sample Center is a robust Autosampler for high sample throughput, optimized for the challenges of the modern analytic laboratory. It works in accordance with the x-y-z principle and a syringe control with high resolution for precise sample injection. A double needle system enables the penetration of vessel caps and septa. A PEEK injection valve completes the automation system of the 889 IC Sample Center. It is the appropriate automation solution when you have only a very small amount of sample.</p>
1 pcs	6.1608.110	Wash bottle 250 mL to IC Sample Center Bottle for washing solution in the IC Sample Center.	
2 pcs	6.1807.010	Y-connector for tubing 6-9 mm i.d. Connector for waste tubing	
2 pcs	6.1816.070	Silicon tubing 8 mm i.d. /1 m For the Sample Center	

1 pcs **6.1831.010** **PEEK capillary 0.25 mm i.d. / 3 m**
For all IC components.



2 pcs **6.2041.200** **Sample rack 48 x 11 mm to IC Sample Center**
Sample rack for the IC Sample Center with 48 sample positions. For sample vials 6.2743.1xx.



1 pcs **6.2151.020** **Cable USB A - USB B 1.8 m**
USB connecting cable



1 pcs **6.2743.107** **PP Vials 0.7 mL to IC Sample Center, 100 pieces**
100 sample vials.



1 pcs **6.2743.127** **Septum to 6.2743.1XX Si/PTFE**
100 pieces










1 pcs **6.2744.010** **Pressure screw 5x**
With UNF 10/32 connection. For the connection of PEEK capillaries



1 pcs **6.2816.200** **Air needle 80 mm, black**
To IC Sample Center



Optional accessories

Order no.	Description	
6.1448.050	Aluminum septum caps / 1000 pieces Including Al sealing.	
6.1448.067	Screw cap with septum, 100 pcs. Screw cap with septum, suitable for sample glass 8 mL 6.2420.007, 100 pieces	
6.1562.170	Tefzel transfer tubing 1000 µL To IC Sample Center.	
6.1829.060	Capillary set for syringe unit For the IC Sample Center	
6.2041.210	Sample rack, 12 x 22 mm, for IC Sample Center Sample rack for 12 vials.	
6.2419.000	Sample glass, 6 mL, 1,000 pieces Sample glass, 6 mL, suitable for 61448050 septum, 1,000 pieces	
6.2420.007	Sample glass with thread, 8 mL, 100 pieces Sample glass with thread, 8 mL, suitable for 6.1448.067 screw cap, 100 pieces	

6.2621.110 Crimping tongs

Crimping tongs for sealing 6 mL sample vials (6.2419.000, 6.2419.007).



6.2743.100 PP Vials 0.7 mL to IC Sample Center, 1000 pieces

1000 sample vials.



6.2743.110 PP Vials 0.3 mL to IC Sample Center, 1000 pieces

1000 sample vials.



6.2743.117 PP Vials 0.3 mL to IC Sample Center, 100 pieces

100 sample vials.



6.2743.120 Septum to 6.2743.1XX Si/PTFE

1000 pieces



6.2816.210 Needle with capillary, 15 μ L

To IC Sample Center



6.5333.000 IQ/OQ Kit for IC

The IQ/OQ Kit contains all parts and standard solutions required for IQ/OQ in ion chromatography.



6.7303.100 Snap-ring bottle ND11, 1.5 mL

For liquid holder with MMS 5000, 100 pcs.



6.7303.110 Septum cap ND11

For liquid holder with MMS 5000, 100 pcs.

