



## Accessories for 2.800.0010

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery 2.800.0010

Qty.	Order no.	Description	
1 pcs	1.800.0030	<b>800 Dosino (cable length 1.5 m)</b>	
		Drive with write/read hardware for intelligent dosing units. With permanently attached cable (length 1.5 m).	
1 pcs	6.1546.030	<b>Extractor rod</b>	
		For PTFE piston of dosing unit	

## Optional accessories

Order no.	Description	
<b>6.1562.130</b>	<b>Transfer Tubing with holder 2 x M6 / 10 mL to Sample Processors</b>	
	Transfer tubing for dilution and preconcentration using a Dosino and Sample Processor.	
<b>6.2151.010</b>	<b>Cable MSB plug / socket</b>	
	Extension cable for MSB connections and for Touch Control - Titrande, Mini-DIN socket - mini-DIN plug	
<b>6.3032.120</b>	<b>Dosing Unit 2 mL</b>	
	Dosing unit with integrated data chip with 2 mL glass cylinder and light protection, mountable on reagent bottle with ISO/DIN GL45 glass thread. FEP tubing connection, anti-diffusion buret tip.	
<b>6.3032.150</b>	<b>Dosing Unit 5 mL</b>	
	Dosing unit with integrated data chip with 5 mL glass cylinder and light protection, mountable on reagent bottle with ISO/DIN GL45 thread. FEP tubing connection, anti-diffusion buret tip.	
<b>6.3032.210</b>	<b>Dosing Unit 10 mL</b>	
	Dosing unit with integrated data chip with 10 mL glass cylinder and light protection, mountable on reagent bottle with ISO/DIN GL45 thread. FEP tubing connection, anti-diffusion buret tip.	
<b>6.3032.220</b>	<b>Dosing Unit 20 mL</b>	
	Dosing unit with integrated data chip with 20 mL glass cylinder and light protection, mountable on reagent bottle with ISO/DIN GL45 thread. FEP tubing connection, anti-diffusion buret tip.	
<b>6.3032.250</b>	<b>Dosing Unit 50 mL</b>	
	Dosing unit with integrated data chip with 50 mL glass cylinder and light protection, mountable on reagent bottle with ISO/DIN GL45 thread. FEP tubing connection, anti-diffusion buret tip.	