







## Accessories for 2.944.0010

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery 2.944.0010

Qty.	Order no.	Description	
1 pcs	1.944.0010	<b>944 Professional UV/VIS Detector Vario</b>	
		The 944 Professional UV/VIS Detector Vario is the intelligent UV/VIS detector from Metrohm. It permits a secure and reliable quantification of substances active in the ultraviolet or visible range. Detection takes place via a Diode Array.	
2 pcs	6.1831.100	<b>PEEK capillary 0.25 mm i.d., 1 m</b>	
1 pcs	6.2151.020	<b>Cable USB A - USB B 1.8 m</b> USB connecting cable	
1 pcs	6.2744.070	<b>Pressure screw short</b> Short version. With UNF 10/32 connection. 5 pieces. For the connection of PEEK capillaries	
1 pcs	6.2839.130	<b>Flow cell 10 mm to UV/VIS instruments</b> PEEK flow cell	

## Optional accessories

Order no.	Description	
<b>2.943.0110</b>	<b>943 Professional Reactor Vario</b>	
	The 943 Professional Reactor Vario is a heatable reactor, especially developed for pre- and post-column derivatization. It captivates by its robustness, the fast heating rate and the possibility of carrying out reactions at up to 150 °C.	
<b>2.943.0210</b>	<b>943 Professional Thermostat Vario</b>	
	The 943 Professional Thermostat Vario is a stand-alone column oven for up to two columns with a maximum length of 150 mm. It captivates by its fast heating rate and the possibility of operating two separation columns together at up to 80°C.	
<b>6.2061.100</b>	<b>Bottle holder for Professional IC instruments</b>	
	Instrument cover and bottle holder for Professional IC instruments.	
<b>6.2061.110</b>	<b>Tray with sensor for Professional IC Instruments</b>	
<b>6.2061.120</b>	<b>System Connector</b>	
	Connects an IC module with the footprint of a Professional IC with a Compact IC	
<b>6.2764.000</b>	<b>Replacement seals to 62839130</b>	
	Set of seals for Flow cell 10 mm to UV/VIS instruments	

---

**6.2804.060 Deuterium lamp (UV) for UV/VIS instruments**

Replacement bulb for the UV range



---

**6.2804.100 Halogen lamp (Vis) for Professional UV/VIS Detector**

Spare lamp for the visible range



---

**6.6059.322 MagIC Net 3.2 Professional CD: 1 license**

Professional PC program for controlling all of the intelligent Professional IC systems, Compact IC instruments and their peripherals, such as various Autosamplers, 800 Dosinos, 771 Compact Interface, etc. The software permits the inspection, data acquisition, data evaluation, data monitoring and report generation of ion chromatographic analyses.



Graphical user interface for routine operations, extensive database programs, method development, configuration and manual system control; very flexible user administration, efficient database operations, extensive data export functions, individually configurable report generator, control and monitoring of all system components and the chromatography results.

MagIC Net Professional complies fully with FDA Regulation 21 CFR Part 11 as well as GLP.

MagIC Net is available in 16 dialog languages.

1 license.

---

**6.6059.323 MagIC Net 3.2 Multi CD: 3 licenses**

Professional PC program for controlling all of the intelligent Professional IC systems, Compact IC instruments and their peripherals, such as various Autosamplers, 800 Dosinos, 771 Compact Interface, etc. The software permits the inspection, data acquisition, data evaluation, data monitoring and report generation of ion chromatographic analyses.



Graphical user interface for routine operations, extensive database programs, method development, configuration and manual system control; very flexible user administration, efficient database operations, extensive data export functions, individually configurable report generator, control and monitoring of all system components and the chromatography results.

MagIC Net Multi complies fully with FDA Regulation 21 CFR Part 11 as well as with GLP.

MagIC Net is available in 16 dialog languages.

Client/server version with 3 licenses.

---