




Accessories for 2.945.0010

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.945.0010

Qty.	Order no.	Description	
1 pcs	1.945.0010	945 Professional Detector Vario – Conductivity	
		Intelligent stand-alone detector equipped with a high-performance IC Conductivity Detector. For use with intelligent IC instruments or as independent conductivity detector.	
1 pcs	6.1831.030	PEEK capillary 0.75 mm i.d. / 3 m	
		For all IC components.	
1 pcs	6.2151.020	Cable USB A - USB B 1.8 m	
		USB connecting cable	
1 pcs	6.2744.014	Pressure screw 2x	
		With UNF 10/32 connection. For the connection of PEEK capillaries	