






## Accessories for 2.942.0070

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery 2.942.0070

Qty.	Order no.	Description	
1 pcs	1.942.0070	<b>942 Extension Module Vario LQH</b>	
		The 942 Extension Module Vario LQH is an extension module for the Professional IC Vario line of instruments. It expands the possibilities of <b>Metrohm Inline Sample Preparation "MISP" and Liquid Handling and opens a Professional IC system for Online Monitoring.</b>	
		Typical applications:	
		<ul style="list-style-type: none"> <li>• <b>Extended Inline Sample Preparation</b></li> <li>• <b>Online monitoring</b></li> <li>• <b>Dilution without sample changers</b></li> <li>• <b>Fully automatic sample derivation reactions</b></li> </ul>	
1 pcs	6.1562.160	<b>Transfer tubing 15 mL, 2 x M6</b>	
		For liquid handling in ion chromatography. To avoid contamination of the Dosino.	
5 pcs	6.1608.100	<b>PE bottle 50 mL</b>	
		Bottle for Extension Module Liquid Handling.	
1 pcs	6.1803.020	<b>PTFE capillary 0.97 mm i.d. / 5 m</b>	
		For all IC instruments	
1 pcs	6.1825.230	<b>PEEK sample loop 10 µL</b>	
		For injection valve, with 2 PEEK pressure screws	

---

**2 pcs    6.1826.340    Pump tubing LFL (black/black), 3 stoppers**

For the sample solution in inline dialysis.



---

**1 pcs    6.1831.010    PEEK capillary 0.25 mm i.d. / 3 m**

For all IC components.



---

**1 pcs    6.2057.200    Bottle holder to Extension Module - Liquid Handling**



---

**1 pcs    6.2156.060    Cable Extension Module - Professional IC, 40 cm**

Cable connecting the extension module with a Professional IC instrument



---

**1 pcs    6.2251.000    Colored sleeves for capillaries**

Colored pieces of heat shrink tubing for capillary coding. Three pieces each of five different colors.



---

**1 pcs    6.2744.010    Pressure screw 5x**

With UNF 10/32 connection. For the connection of PEEK capillaries



---

**1 pcs    6.2744.034    Coupling olive/UNF 10/32, 2x**

Connection between pressure screw and pump tubing. 2 pieces. For IC instruments with peristaltic pump.



---

**2 pcs    6.2744.070    Pressure screw short**

Short version. With UNF 10/32 connection. 5 pieces. For the connection of PEEK capillaries



---

**2 pcs    6.2744.080    M6 thread / UNF 10/32 adapter**

With M6 outer thread. For IC instruments.



---

**6 pcs    6.2744.090    Pressure screw long**

Long version. With UNF 10/32 connection. 2 pieces. For the connection of PEEK capillaries (MCS and Sample degasser)



---

**2 pcs    6.2744.160    Nozzle/UNF 10/32 coupling with security device**

Security device for connecting pump tubing to a nozzle.





---

**1 pcs**      **6.2744.200**      **Adapter UNF 10/32 outer / M6 inner**  
Connecting 1/16 in. PEEK capillaries and transfer tubing (M6)



## Optional accessories

Order no.	Description	
<b>6.2103.170</b>	<b>Adapter to Leak Sensor at Professional IC</b> Adapter to connect two leak sensors to a Professional IC instrument	
<b>6.2156.070</b>	<b>Cable Extension Module - Professional IC, 1 m</b> Cable connecting the extension module with the Professional IC instrument	
<b>6.5332.000</b>	<b>IC tubing kit for second extension module</b>	