





Accessories for 2.1001.2220

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.1001.2220

Qty.	Order no.	Description	
1 pcs	1.1001.0020	OMNIS Titrator with stirrer	
		OMNIS Titrator with magnetic stirrer, without function license	
1 pcs	6.00201.300	dEcotrode plus	
		Digital, combined pH electrode for OMNIS. This electrode is suitable for aqueous acid/base titrations. The fixed ground-joint diaphragm is insensitive to contamination. Reference electrolyte: c(KCl) = 3 mol/L, storage in storage solution.	
1 pcs	6.01600.010	OMNIS Liquid Adapter	
		Adapter for connecting reagent bottles with 3S technology to an OMNIS Titrator, Titration Module or Dosing Module.	
1 pcs	6.01601.000	Bottle cap, multi-use	
		Reusable bottle cap with 3S technology for connecting reagent bottles to an OMNIS Titrator, Titration Module or Dosing Module.	

1 pcs	6.02004.010	Bottle holder		
			Bottle holder for OMNIS Titrator, Titration Module and Dosing Module.	
1 pcs	6.02005.000	OMNIS electrode holder		
			Electrode holder for one OMNIS titration system.	
1 pcs	6.02005.010	Electrode Parking Station		
			Electrode Parking Station for electrode storage on an OMNIS titration system.	
1 pcs	6.02100.010	Digital measuring module		
			Digital measurement channel for one OMNIS Titrator or Titration Module for the connection of digital "dTrodes" electrodes.	
1 pcs	6.02104.300	Electrode cable plug-in head Q / plug P, 0.55 m		
			Cable for connecting a dTrobe (plug-in head Q) to the digital measuring module from OMNIS.	
1 pcs	6.02107.000	Network cable, 1 m		
			Network cable (1 m) for connecting one OMNIS Sample Robot or OMNIS Titrator to an Ethernet Router, Ethernet Switch or local network	

1 pcs 6.03001.220 OMNIS 20 mL cylinder unit

Intelligent 20 mL cylinder unit for an OMNIS Titrator, Titration Module or Dosing Module. Includes dosing tubing and antidiffusion tip.



1 pcs 6.06000.000 OMNIS Software storage device

License-based PC program for controlling complex OMNIS analysis systems.

- Sample-oriented operating procedures with numerous templates
- Graphical method editor with numerous templates
- Clearly organized and professional sample lists
- Management of instruments, work systems, sensors, solutions
- Dialog languages: German, English, Chinese



OMNIS
A WHOLE NEW LEVEL OF PERFORMANCE

1 pcs 6.06001.910 Function license Advanced

"Advanced" function license for one OMNIS Titrator

- Liquid Handling (LQH)
- MEAS U/T/pH
- Calibration (Cal)
- Endpoint titration without conditioning (SET)
- Monotonic equivalence point titration (MET)
- Dynamic equivalence point titration (DET)



OMNIS
A WHOLE NEW LEVEL OF PERFORMANCE

1 pcs 6.06003.010 OMNIS Stand-Alone license

Enables stand-alone operation of the OMNIS software on a Windows™ computer.

Features:

- The license already includes one OMNIS instrument license.
- Must be activated via the Metrohm licensing portal.
- Not transferable to another computer.



OMNIS
A WHOLE NEW LEVEL OF PERFORMANCE

1 pcs 6.1546.030 Extractor rod

For PTFE piston of dosing unit



1 pcs **6.1608.020** **Glass bottle / 1000 mL / GL 45**

Bottle for auxiliary solutions.



1 pcs **6.1903.020** **Stirring bar / 16 mm**

Stirring Bar with magnetic core, PTFE coated, length 16 mm



1 pcs **6.1903.030** **Stirring bar / 25 mm**

Stirring bar with magnetic core, PTFE covering, length 25 mm.



1 pcs **6.2013.010** **Clamping ring**

For support rods with a diameter of 10 mm.



1 pcs **6.2016.050** **Support rod / 300 mm**



1 pcs **6.2307.230** **Buffer solutions pH 4, 7 and 9**

Mixed buffer solutions pH 4.00/7.00/9.00 (25 °C) in single use sachets, colourless, box of 3 x 10 x 30 mL



1 pcs **6.2325.000** **pHit kit**




Maintenance kit for pH electrodes

The kit contains:

- 50 mL cleaning solution
- 50 mL 3M KCl solution
- 50 mL storage solution
- 2 Storage vessels
- Instructions for use



Optional accessories

Order no.	Description	
6.05609.000	KFT equipment for OMNIS Titrator	
	Complete accessory set for volumetric Karl Fischer titration.	
6.05609.010	KFT equipment for OMNIS Titrator with Polytron	
	Complete accessory set for volumetric Karl Fischer titration with Polytron.	
6.06001.160	Function license KFT with conditioning	
	Function license for volumetric Karl Fischer titration with conditioning for an OMNIS titration system.	