





## Accessories for 2.1003.0010

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery 2.1003.0010

Qty.	Order no.	Description	
1 pcs	1.1003.0010	<b>OMNIS Dosing Module</b>	
		Dosing module for extending an OMNIS titration system to include an additional buret.	
1 pcs	6.01600.010	<b>OMNIS Liquid Adapter</b>	
		Adapter for connecting reagent bottles with 3S technology to an OMNIS Titrator, Titration Module or Dosing Module.	
1 pcs	6.01601.000	<b>Bottle cap, multi-use</b>	
		Reusable bottle cap with 3S technology for connecting reagent bottles to an OMNIS Titrator, Titration Module or Dosing Module.	
1 pcs	6.02004.010	<b>Bottle holder</b>	
		Bottle holder for OMNIS Titrator, Titration Module and Dosing Module.	

---

**1 pcs**      **6.1546.030**                      **Extractor rod**  
For PTFE piston of dosing unit







---

**1 pcs**      **6.1608.020**                      **Glass bottle / 1000 mL / GL 45**  
Bottle for auxiliary solutions.



## Optional accessories

Order no.	Description	
<b>2.1002.0010</b>	<b>OMNIS Titration Module without stirrer</b>	
	<p>Titration module for connection to an OMNIS Titrator for extending the system to include an additional buret for titration/dosing and up to two measuring modules (analog or digital). Can be supplemented with one magnetic stirrer or rod stirrer for use as separate titration stand. Freely selectable cylinder unit with 5, 10, 20 or 50 mL.</p>	
<b>2.1002.0110</b>	<b>OMNIS Titration Module with magnetic stirrer</b>	
	<p>Titration module with magnetic stirrer for connection to an OMNIS Titrator for extending the system to include an additional titration stand, a buret for titration/dosing and up to two measuring modules (analog or digital). Freely selectable cylinder unit with 5, 10, 20 or 50 mL.</p>	
<b>2.1003.0010</b>	<b>OMNIS Dosing Module without stirrer</b>	
	<p>Dosing module for connection to an OMNIS Titrator for extending the system to include an additional buret for titration/dosing. Can be supplemented with one magnetic stirrer or rod stirrer for use as separate titration stand. Freely selectable cylinder unit with 5, 10, 20 or 50 mL.</p>	
<b>2.1003.0110</b>	<b>OMNIS Dosing Module with magnetic stirrer</b>	
	<p>Dosing module with magnetic stirrer for connection to an OMNIS Titrator for extending the system to include an additional titration stand and a buret for titration/dosing. Freely selectable cylinder unit with 5, 10, 20 or 50 mL.</p>	
<b>6.01503.120</b>	<b>OMNIS 2 mL cylinder unit without accessories</b>	
	<p>Intelligent 2 mL cylinder unit for an OMNIS Titrator, Titration Module or Dosing Module. Without dosing tubing and antidiffusion tip.</p>	

---

**6.01503.150      OMNIS 5 mL cylinder unit without accessories**

Intelligent 5 mL cylinder unit for an OMNIS Titrator, Titration Module or Dosing Module. Without dosing tubing and antidiffusion tip.



---

**6.01503.210      OMNIS 10 mL cylinder unit without accessories**

Intelligent 10 mL cylinder unit for an OMNIS Titrator, Titration Module or Dosing Module. Without dosing tubing and antidiffusion tip.



---

**6.01503.220      OMNIS 20 mL cylinder unit without accessories**

Intelligent 20 mL cylinder unit for an OMNIS Titrator, Titration Module or Dosing Module. Without dosing tubing and antidiffusion tip.



---

**6.01503.250      OMNIS 50 mL cylinder unit without accessories**

Intelligent 50 mL cylinder unit for an OMNIS Titrator, Titration Module or Dosing Module. Without dosing tubing and antidiffusion tip.



---

**6.02102.010      Cable MDL PL/SO 0.5 m**

Connecting cable for OMNIS instruments with MDL connector (0.5 m)



---

**6.02102.020**      **Cable MDL PL/SO 1 m**

Connecting cable for OMNIS instruments with MDL connector (1 m)



---

**6.02102.030**      **Cable MDL PL/SO 2 m**

Connecting cable for OMNIS instruments with MDL connector (2 m)



---

**6.02102.040**      **Cable MDL PL/SO 5 m**

Connecting cable for OMNIS instruments with MDL connector (5 m)



---

**6.02102.050**      **Cable MDL PL/SO 10 m**

Connecting cable for OMNIS instruments with MDL connector (10 m)



---

**6.02106.000**      **Liquid Adapter extension cable, 0.55 m**

Extension cable (0.55 m in length) for connecting an OMNIS Liquid Adapter to an OMNIS Titrator, Titration Module or Dosing Module.



---

**6.02106.010**      **Liquid Adapter extension cable, 2 m**

Extension cable (2 m in length) for connecting an OMNIS Liquid Adapter to an OMNIS Titrator, Titration Module or Dosing Module.



---

**6.03001.120**      **OMNIS 2 mL cylinder unit**

Intelligent 2 mL cylinder unit for an OMNIS Titrator, Titration Module or Dosing Module. Includes dosing tubing and antidiffusion tip.



---

**6.03001.150      OMNIS 5 mL cylinder unit**

Intelligent 5 mL cylinder unit for an OMNIS Titrator, Titration Module or Dosing Module. Includes dosing tubing and antidiffusion tip.



---

**6.03001.210      OMNIS 10 mL cylinder unit**

Intelligent 10 mL cylinder unit for an OMNIS Titrator, Titration Module or Dosing Module. Includes dosing tubing and antidiffusion tip.



---

**6.03001.220      OMNIS 20 mL cylinder unit**

Intelligent 20 mL cylinder unit for an OMNIS Titrator, Titration Module or Dosing Module. Includes dosing tubing and antidiffusion tip.



---

**6.03001.250      OMNIS 50 mL cylinder unit**

Intelligent 50 mL cylinder unit for an OMNIS Titrator, Titration Module or Dosing Module. Includes dosing tubing and antidiffusion tip.



---

**6.05609.000      KFT equipment for OMNIS Titrator**

Complete accessory set for volumetric Karl Fischer titration.



---

**6.05609.010      KFT equipment for OMNIS Titrator with Polytron**

Complete accessory set for volumetric Karl Fischer titration with Polytron.

