



## Accessories for 2.800.0020

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery 2.800.0020

Qty.	Order no.	Description	
1 pcs	1.800.0040	<b>800 Dosino (cable length 0.65 m)</b>	
		Drive with write/read hardware for intelligent dosing units. With permanently attached cable (length 0.65 m).	
1 pcs	6.1546.030	<b>Extractor rod</b>	
		For PTFE piston of dosing unit	