

6.1115.000

6.1115.000

1	6.1903.030	Stirring bar / 25 mm
----------	-------------------	-----------------------------

Stirring bar with magnetic core, PTFE covering, length 25 mm.



1	6.2008.040	Storage vessel
----------	-------------------	-----------------------

Together with 6.2008.050 storage vessel holder provides a support for the electrode on 6.3032.XXX dosing units.



6.02109.000**OMNIS adapter cable**

This adapter cable comprising socket F to plug P is suitable for use as a fixed cable (e.g. Optrode or Solvotrode easyClean) connecting a sensor without temperature sensor to the OMNIS analog measuring module.



6.2166.000**USB power supply unit 5V 1A R**

Power supply unit for the power supply of the 6.1115.100 Optrode if there is no USB connector available for this purpose on a Metrohm device (e.g. with older devices).

