

6.3032.150

6.3032.150

-
- | | | |
|----------|-------------------|---------------------------------|
| 1 | 6.1446.030 | Ball stopper / B-14/(15) |
|----------|-------------------|---------------------------------|
- For the movable holding of buret tips in titration vessels



-
- | | | |
|----------|-------------------|------------------------------|
| 1 | 6.1446.040 | Threaded stopper / M6 |
|----------|-------------------|------------------------------|



-
- | | | |
|----------|-------------------|------------------------|
| 1 | 6.1543.060 | Tip / M6 thread |
|----------|-------------------|------------------------|



Tip with M6 thread. Together with the 6.1805.160 capillary and the 6.1446.030 link stopper, it makes up the entire 6.1537.010 buret tip.



-
- | | | |
|----------|-------------------|----------------------------------|
| 1 | 6.1543.200 | Titration tip / M6 thread |
|----------|-------------------|----------------------------------|

Titration tip with antidiffusion valve.



1	6.1619.000	Adsorber tube for Dosing Unit	
1	6.1805.100	FEP tubing / M6 / 40 cm With light and kink protection	
1	6.1829.010	FEP aspiration tube M6, 0.25 m With M 6 thread. For Dosing Unit	
1	6.2008.030	Buret tip storage vessel for Dosing Unit Together with 6.2008.050 buret tip holder provides a fastening for the buret tips on 6.3032.XXX Dosing units.	

-
- 1** **6.2008.050** **Storage vessel holder for Dosing Unit**
Together with 6.2008.030 buret tip storage vessel provides a fastening for the buret tip on Dosing units.



-
- 1** **6.2244.020** **Cards for intelligent Exchange Unit**
Set of 10 pieces, different colors.



-
- 1** **6.2739.000** **Wrench**
For tightening connectors



-
- 0** **6.2803.000** **Grease (1 oz.)**
Special quality (without silicone). For grease-free ground joint connections see 6.2713.XXX



-
- 1** **6.2803.010** **Grease (2 g)**
Special quality (without silicone). For grease-free ground joint connections see 6.2713.XXX



6.1608.020 Glass bottle / 1000 mL / GL 45

Bottle for auxiliary solutions.



6.1608.023 Amber glass bottle / 1000 mL / GL 45

For Exchange Units. Bottle for auxiliary solutions.



6.1608.030 Round glass bottle / 1000 mL / GL 45

6.1608.040 PE bottle / 1000 mL / GL 45

For Exchange Units. Bottle for auxiliary solutions.



6.1608.050 Drying bottle / 100 mL / GL 45

6.1608.060**Bottle / 100 mL**

6.1618.000**Thread adapter 32 mm to GL 45**

Adapter for bottles with 32 mm thread (Fluka, Riedel de-Haen 500 mL...) to GL45




6.1618.010**Thread adapter 28 mm to GL 45**

Adapter for bottles with 28 mm thread (Fisher...) to GL45



6.1618.020**Thread adapter S 40 to GL 45**

Adapter for bottles with S40 thread (Merck...) to GL45.



6.1618.050**Thread adapter 40 mm to GL 45**

Adapter for bottles with 40 mm thread to GL45.



6.1618.060**Thread adapter GL 45 to GL 45**

Thread adapter GL45 to GL45; e.g. for reagent bottles with extended flange

