

2.945.0010

2.945.0010

1	1.945.0010	945 Professional Detector Vario – Conductivity
----------	-------------------	---

Intelligent stand-alone detector equipped with a high-performance IC Conductivity Detector. For use with intelligent IC instruments or as independent conductivity detector.



1	6.1831.030	PEEK capillary 0.75 mm i.d. / 3 m
----------	-------------------	--

For all IC components.



1	6.2151.020	Cable USB A - USB B 1.8 m
----------	-------------------	----------------------------------

USB connecting cable



1	6.2744.014	Pressure screw 2x
----------	-------------------	--------------------------

With UNF 10/32 connection. For the connection of PEEK capillaries



6.6059.322 MagIC Net 3.2 Professional CD: 1 license



Professional PC program for controlling all of the intelligent Professional IC systems, Compact IC instruments and their peripherals, such as various Autosamplers, 800 Dosinos, 771 Compact Interface, etc. The software permits the inspection, data acquisition, data evaluation, data monitoring and report generation of ion chromatographic analyses.

Graphical user interface for routine operations, extensive database programs, method development, configuration and manual system control; very flexible user administration, efficient database operations, extensive data export functions, individually configurable report generator, control and monitoring of all system components and the chromatography results.

MagIC Net Professional complies fully with FDA Regulation 21 CFR Part 11 as well as GLP.

MagIC Net is available in 16 dialog languages.

1 license.

6.6059.323 MagIC Net 3.2 Multi CD: 3 licenses



Professional PC program for controlling all of the intelligent Professional IC systems, Compact IC instruments and their peripherals, such as various Autosamplers, 800 Dosinos, 771 Compact Interface, etc. The software permits the inspection, data acquisition, data evaluation, data monitoring and report generation of ion chromatographic analyses.

Graphical user interface for routine operations, extensive database programs, method development, configuration and manual system control; very flexible user administration, efficient database operations, extensive data export functions, individually configurable report generator, control and monitoring of all system components and the chromatography results.

MagIC Net Multi complies fully with FDA Regulation 21 CFR Part 11 as well as with GLP.

MagIC Net is available in 16 dialog languages.

Client/server version with 3 licenses.
