

2.800.0010

2.800.0010

1	1.800.0030	800 Dosino (cable length 1.5 m)
----------	-------------------	--

Drive with write/read hardware for intelligent dosing units. With permanently attached cable (length 1.5 m).



1	6.1546.030	Extractor rod
----------	-------------------	----------------------

For PTFE piston of dosing unit



6.1562.130 Transfer Tubing with holder 2 x M6 / 10 mL to Sample Processors

Transfer tubing for dilution and preconcentration using a Dosino and Sample Processor.



6.2151.010 Cable MSB plug / socket

Extension cable for MSB connections and for Touch Control - Titrande, Mini-DIN socket - mini-DIN plug



6.3032.120 Dosing Unit 2 mL

Dosing unit with integrated data chip with 2 mL glass cylinder and light protection, mountable on reagent bottle with ISO/DIN GL45 glass thread. FEP tubing connection, anti-diffusion buret tip.



6.3032.150 Dosing Unit 5 mL

Dosing unit with integrated data chip with 5 mL glass cylinder and light protection, mountable on reagent bottle with ISO/DIN GL45 thread. FEP tubing connection, anti-diffusion buret tip.



6.3032.210 Dosing Unit 10 mL

Dosing unit with integrated data chip with 10 mL glass cylinder and light protection, mountable on reagent bottle with ISO/DIN GL45 thread. FEP tubing connection, anti-diffusion buret tip.



6.3032.220 Dosing Unit 20 mL

Dosing unit with integrated data chip with 20 mL glass cylinder and light protection, mountable on reagent bottle with ISO/DIN GL45 thread. FEP tubing connection, anti-diffusion buret tip.



6.3032.250 Dosing Unit 50 mL

Dosing unit with integrated data chip with 50 mL glass cylinder and light protection, mountable on reagent bottle with ISO/DIN GL45 thread. FEP tubing connection, anti-diffusion buret tip.

