






Accessories for 2.915.0110

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.915.0110

Qty.	Order no.	Description	
1 pcs	1.800.0040	800 Dosino (cable length 0.65 m)	
		Drive with write/read hardware for intelligent dosing units. With permanently attached cable (length 0.65 m).	
1 pcs	1.915.0010	915 KF Ti-Touch	
1 pcs	6.0338.100	Double Pt-wire electrode for volumetry	
		Indicator electrode used for volumetric Karl Fischer titration.	
1 pcs	6.1244.040	Set of seals	
		Set of seals for 6.1414.030. The set is comprised of one seal and various O-rings.	
1 pcs	6.1403.040	KF Adsorber tube	
		Used with KF equipment. Including cover and O-ring.	

1 pcs **6.1414.030** **KF titration vessel lid with 5 openings**
For volumetric KF titrations.



1 pcs **6.1415.220** **Titration vessel / 20-90 mL**



1 pcs **6.1415.250** **Titration vessel / 50-150 mL**



2 pcs **6.1446.040** **Threaded stopper / M6**



1 pcs **6.1446.090** **Stopper / B-14/15 / M8**



2 pcs **6.1448.010** **Septum 12 mm, 5 pieces**
Set of 5 items.



1 pcs **6.1543.110** **Dosing tip / M8 thread**
Transfer tip with M8 thread



1 pcs **6.1546.030** **Extractor rod**
For PTFE piston of dosing unit








2 pcs **6.1602.105** **Bottle attachment / GL 45 / green**
For direct dosing from reagent bottles with GL 45 thread.



1 pcs **6.1608.023** **Amber glass bottle / 1000 mL / GL 45**
For Exchange Units. Bottle for auxiliary solutions.



1 pcs	6.1608.030	Round glass bottle / 1000 mL / GL 45	
1 pcs	6.1609.010	Adsorber tube with tubing nipple Lid with tubing nipple.	
1 pcs	6.1623.000 For 803 Ti Stand	Overflow protection for Ti Stand	
1 pcs	6.1801.120	PVC tubing / 4 mm / 6 mm / 2 m	
2 pcs	6.1805.200	PTFE tubing / M8 / 0.5 m	

1 pcs **6.1808.050** **Tubing adapter olive / M8 outer**
1 M8 outer thread and 1 tubing olive. E.g. for thermostat jacket of exchange units and stability measuring instruments.



1 pcs **6.1819.030** **Tube**
For M8 threaded holes. For 6.1602.1XX bottle attachments and 6.1608.0XX bottles



1 pcs **6.1819.050** **Tube**
For M8 threaded holes of 6.1602.10X Bottle attachments with 6.1608.030 Aspiration bottle. For 703 and 803 Ti Stand, 681 Pump Unit and 679 Rancimat



1 pcs **6.1821.070** **Aspiration tip with M10 screw nipple**
For aspirating reagent in combination with 803 Ti Stand and volumetric Karl Fischer cell. The aspiration tip fits tubing with an M8 thread. The attached tubing nipple has an M10 thread.



2 pcs **6.1903.020** **Stirring bar / 16 mm**
Stirring Bar with magnetic core, PTFE coated, length 16 mm



2 pcs **6.1903.030** **Stirring bar / 25 mm**
Stirring bar with magnetic core, PTFE covering, length 25 mm.



1 pcs **6.2013.010** **Clamping ring**
For support rods with a diameter of 10 mm.



1 pcs **6.2016.050** **Support rod / 300 mm**



1 pcs **6.2023.020** **Clip for SGJ 14/15**
Clip for SGJ 14/16



2 pcs **6.2043.005** **Holding clip for bottles**

Holding spring for reagent bottles in exchange units.



1 pcs **6.2104.020** **Electrode cable / 1 m / F**

For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).



1 pcs **6.2164.010** **Power supply unit 100 - 240 V/24 V DC**

Power supply unit for 915 KF Ti-Touch and 916 Ti-Touch



1 pcs **6.2412.000** **Weighing boats**

Weighing boat made of glass for the 6.1414.030 titration vessel lid or the 6.1455.31X, 6.1464.32X, 6.1465.320 titration vessel. For KF water determinations. Including protective tube.



1 pcs **6.2621.030** **Hexagon key 4 mm**



2 pcs **6.2730.010** **Screw nipple**

Made of PP. For 6.1448.010 Septum or for mounting 6.1403.040 Drying tube or 6.0338.100 Electrode in 6.1414.030 and 6.1414.040 KF Titration vessel lids. Used for KF titrations.



1 pcs 6.2730.020 Septum stopper

With E.301.0041 O-ring. For 6.1448.010 Septum and 6.1414.030 Titration vessel lid. Used for KF applications.



3 pcs 6.2730.030 Stopper

Complete with nipple and O-ring. For 6.1414.030 titration vessel lids. For KF instruments, VA measurement stands, VA Computrace and IC Dilution.



1 pcs 6.2739.000 Wrench

For tightening connectors



1 pcs 6.2811.000 Molecular sieve

Molecular sieve. Bottle containing 250 g. Pore size: 0.3 nm. Without moisture indicator. For Rancimats and Karl Fischer instruments.



1 pcs 6.3032.210 Dosing Unit 10 mL

Dosing unit with integrated data chip with 10 mL glass cylinder and light protection, mountable on reagent bottle with ISO/DIN GL45 thread. FEP tubing connection, anti-diffusion buret tip.



1 pcs **6.6066.000** **USB flash drive for Ti-Touch**
USB flash drive for 915 KF Ti-Touch and 916 Ti-Touch.



Optional accessories

Order no.	Description
-----------	-------------

2.136.0100	Polytron PT 1300 D
-------------------	---------------------------

Polytron PT 1300 D - Metrohm version

Homogenizer that can be controlled directly by OMNIS Software, tiamo™ or Touch Control.

The Polytron PT 1300 D is comprised of a control unit and a drive. The coupling system on the drive permits simple and rapid replacement of the aggregates without additional tools.

Solid samples can be pulverized without difficulty. However the instrument is also highly suitable for ensuring thorough mixing in viscous samples.



2.142.0100	Custom Q3X thermal printer
-------------------	-----------------------------------

Compact printer with USB interface for

- 900 Touch Control
- 915 KF Ti-Touch
- 916 Ti-Touch
- 917 Coulometer
- 877 / 848 Titrino plus
- 865 / 876 Dosimat plus
- 91X Meter (cable 6.2151.140)
- Eco Dosimat / Titrator
- 862 Compact Titrosampler
- 870 KF Titrino plus
- 899 Coulometer

Paper width 60 mm (40 characters). Including 6.2151.120 USB cable.



2.800.0010	800 Dosino
-------------------	-------------------

Drive with write/read hardware for intelligent Dosing Units. With fixed cable (length 150 cm).



2.805.0010	805 Dosimat
-------------------	--------------------

Dosing device for the Titrando and Sample Processors with read/write hardware for intelligent Exchange Units. With permanently attached cable. Without Exchange Unit.



2.810.0010 **810 Sample Processor**

Sample Processor for the automated analysis of routine samples using 916 Ti-Touch or 915 KF Ti-Touch.

Sample Processor with one workstation and one built-in membrane pump for the automatic processing of potentiometric as well as Karl Fischer titrations, in series with small to medium quantities. In addition to the built-in pump, an additional one (membrane or peristaltic) and up to three dosing devices for Liquid Handling tasks can be connected.

Because of the multitude of application variants, rack, stirrer, titration head, Swing Head and sample vessels must be tailored to the application and ordered separately.



6.2047.110 **Holder for Polytron**

Holder to fix the Polytron at the stand rod.



6.2723.300 **Protective cover for 900 Touch Control**

Protective cover for 900 Touch Control



6.3026.110 **Exchange Unit 1 mL**

Exchange Unit with integrated data chip with 1 mL glass cylinder and light protection. PCTFE/PTFE flat stopcock, FEP tubing connection, antidiffusion buret tip and standard amber glass reagent bottle.



6.3026.150 **Exchange Unit 5 mL**

Exchange Unit with integrated data chip with 5 mL glass cylinder and light protection. PCTFE/PTFE flat stopcock, FEP tubing connection, anti-diffusion buret tip and standard amber glass reagent bottle.



6.3026.210 Exchange Unit 10 mL

Exchange unit with integrated data chip with 10 mL glass cylinder and light protection. PCTFE/PTFE flat stopcock, FEP tubing connection, anti-diffusion buret tip and standard amber glass reagent bottle.



6.3026.220 Exchange Unit 20 mL

Exchange unit with integrated data chip with 20 mL glass cylinder and light protection. PCTFE/PTFE flat stopcock, FEP tubing connection, anti-diffusion buret tip and standard amber glass reagent bottle.



6.3026.250 Exchange unit 50 mL

Exchange unit with integrated data chip with 50 mL glass cylinder and light protection. PCTFE/PTFE flat stopcock, FEP tubing connection, antidiiffusion tip and standard amber glass reagent bottle.



6.3032.120 Dosing Unit 2 mL

Dosing unit with integrated data chip with 2 mL glass cylinder and light protection, mountable on reagent bottle with ISO/DIN GL45 glass thread. FEP tubing connection, anti-diffusion buret tip.



6.3032.150 Dosing Unit 5 mL

Dosing unit with integrated data chip with 5 mL glass cylinder and light protection, mountable on reagent bottle with ISO/DIN GL45 thread. FEP tubing connection, anti-diffusion buret tip.



6.3032.220**Dosing Unit 20 mL**

Dosing unit with integrated data chip with 20 mL glass cylinder and light protection, mountable on reagent bottle with ISO/DIN GL45 thread. FEP tubing connection, anti-diffusion buret tip.



6.3032.250**Dosing Unit 50 mL**

Dosing unit with integrated data chip with 50 mL glass cylinder and light protection, mountable on reagent bottle with ISO/DIN GL45 thread. FEP tubing connection, anti-diffusion buret tip.

