


Accessories for 2.136.0400




Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.136.0400

Qty.	Order no.	Description
1 pcs	1.136.0400	PILS, Particle Into Liquid Sampler Particle into Liquid Sampler for the collection of ambient aerosols. Coupling with ion chromatographs allows automatic determination of the ionic composition of the aerosols without requiring additional sample preparation.



Optional accessories

Order no.	Description	
2.136.0500	Peristaltic pump, 8-channel	
	8-channel peristaltic pump for PILS (Particle Into Liquid sampler) for transferring any auxiliary solution.	
6.5335.000	IC equipment for Liquid Handling set for PILS	
	Accessory set for PILS (Particle Into Liquid Sampler).	
6.9014.000	Debubler for PILS	
	To remove the air bubbles prior to analysis.	
6.9014.100	Impactor for PILS including O-ring	