




## Accessories for 2.700.0020

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery 2.700.0020

Qty.	Order no.	Description	
1 pcs	1.700.0020	<b>700 Dosino</b> The 700 Dosino is a Dosing Unit drive. With fixed cable and mini-DIN plug.	
1 pcs	6.1546.030	<b>Extractor rod</b> For PTFE piston of dosing unit	
1 pcs	6.2148.040	<b>Adapter for Dosino 700 (series 23 or higher)</b>	

## Optional accessories

Order no.	Description
<b>6.2134.010</b>	<b>Adapter DB9 male / Mini-DIN8 female</b>
<p>DB9 plug and mini-DIN 8 socket. E.g. for the connection of 700.0020 Dosinos with mini-DIN 8 plugs to 726, 796 Titroprocessors.</p>	

