




Accessories for 2.786.0020

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.786.0020

Qty.	Order no.	Description	
1 pcs	1.786.0010	786 Swing Head with transfer head left Additional equipment for Robotic Sample Processors XL for transferring / pipetting from small sample vessels into larger titration vessels on the rack or into an external titration cell; this can be located to the left of the work station.	
1 pcs	6.1462.040	Robotic arm with transfer head for 786 Swing Head, right swinging Transfer robotic arm for liquid handling purposes	
1 pcs	6.1562.100	Pipetting tubing / 10 mL Pipetting tubing for fully automatic transfer of sample solutions.	
1 pcs	6.1823.010	Guide shaft to pipetting tubing For assembly in the transfer robotic arm on the Robotic Sample Processors	

1 pcs **6.2621.120** **Hexagon key 1.5 mm**

