

## Accessories for 2.786.0030

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery 2.786.0030

Qty.	Order no.	Description	
1 pcs	1.786.0010	<b>786 Swing Head with transfer head left</b> Additional equipment for Robotic Sample Processors XL for transferring / pipetting from small sample vessels into larger titration vessels on the rack or into an external titration cell; this can be located to the left of the work station.	
1 pcs	6.1462.050	<b>Robotic arm with titration head for 786 Swing Head, left or right swinging</b> Titration head for Swing Head	
1 pcs	6.1808.040	<b>Thread adapter M6 outer / M8 inner</b> Outer thread M6, inner thread M8.	
1 pcs	6.1909.030	<b>Stirring propeller / 104 mm</b> Propeller stirrer with PP shaft, PVDF tip. Length from the SGJ: 104 mm. For KF applications with 722, 802 Rod Stirrer	

---

**2 pcs**    **6.2042.030**    **Buret tip clip (rod stirrer)**

Clip for proper alignment of buret tips, dosing or aspiration tubings using a rod stirrer in the titration beaker of a sample changer



---

**1 pcs**    **6.2621.120**    **Hexagon key 1.5 mm**

