



## Accessories for 2.786.0040

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery 2.786.0040

Qty.	Order no.	Description	
1 pcs	1.786.0010	<b>786 Swing Head with transfer head left</b> Additional equipment for Robotic Sample Processors XL for transferring / pipetting from small sample vessels into larger titration vessels on the rack or into an external titration cell; this can be located to the left of the work station.	
1 pcs	6.2621.120	<b>Hexagon key 1.5 mm</b>	

## Optional accessories

Order no.	Description	
<b>6.1462.240</b>	<b>Transfer head for 786 Swing Head, bent and right swinging</b>	
	Transfer head, bent and right swinging, for 786 Swing Head. In combination with the 6.1808.220 Adapter the transfer head can be used on multi-row racks to pick up tools with luer connection.	
<b>6.1462.250</b>	<b>Swing arm for Polytron and 786 Swing Head, left swinging</b>	
	Swing arm for Polytron and 786 Swing Head, left swinging, with one buret tip for the addition of solvent or titrant and three spray nozzles for cleaning. The swing arm allows the usage of the Polytron for sample preparation on multi row racks.	
<b>6.1462.260</b>	<b>Robotic arm with holder for titration head, left swinging, external</b>	
	Titration head holder for 786 Swing Head with possibility of swinging to external positions	
<b>6.2058.040</b>	<b>Swing arm enhancement for Swing Head 786</b>	
	The swing arm enhancement is for stabilizing the swing arm connected to the 786 Swing Head and working with various tools	

---

**6.2058.090      Reinforcement for SwingHead R**

Reinforcement for 786 Swing Head (series 31 and higher) for usage of swing arms with stripping function

