







Accessories for 6.01402.000

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Optional accessories

| Order no. | Description | |
|--------------------|--|---|
| 6.01403.010 | Titration head 4xM10 (P&P) |  |
| | Titration head with 4 M10 openings for use in Pick&Place modules of the OMNIS Sample Robot Pick&Place | |
| 6.01403.020 | Titration head 4xM10 for KFT (P&P) |  |
| | Titration head with 4 M10 openings for Karl Fischer titration in Pick&Place modules of the OMNIS Sample Robot Pick&Place with 4 M10 screw nipples and stoppers. | |
| 6.01403.040 | Titration head 2xNS14 (P&P) |  |
| | Titration head with 2 NS14 openings for use in Pick&Place modules of the OMNIS Sample Robot Pick&Place with 3 retracted dosing tips. Prepared for being equipped with rinsing nozzles. One aspiration tube is already retracted into the titration head. | |
| 6.01404.040 | Beaker adapter for OMNIS 75 mL sample beaker |  |
| | Beaker adapter for centering the 75 mL sample beaker in the slide of the Pick&Place modules | |
| 6.01900.020 | Stirring propeller 13 mm ETFE |  |
| | Stirring propeller 13 mm ETFE with adjustable height for use with Rod Stirrer "Titrator" or Rod Stirrer "Sample Robot" | |
| 6.02041.040 | OMNIS sample rack 25 x 75 mL |  |
| | OMNIS sample rack 25 x 75 mL for use with 75 mL glass sample beaker (P&P) | |

6.02710.000 KF Dis-Cover lid for OMNIS 75 mL sample beaker, 25 pieces

Lid with seals for volumetric Karl Fischer Titration in 6.01402.000 and 6.01402.003 sample beakers in the OMNIS Sample Robot Pick&Place system.



6.02710.030 Dis-Cover lid for OMNIS 75 mL sample beaker, 25 pieces

Dis-Cover lid for 75 mL 6.01402.000 and 6.01402.003 sample beakers in the OMNIS Sample Robot Pick&Place system.

