

## Accessories for 6.01403.000

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery 6.01403.000

Qty.	Order no.	Description	
2 pcs	6.1236.020	<b>Sleeve with SGJ 14/12 mm</b> Sleeve with SGJ 14/12 mm and O-ring.	
5 pcs	6.1446.000	<b>SGJ stopper / B-14/(15)</b>	
3 pcs	6.1446.010	<b>Stopper</b> Used with Sample Changers	
1 pcs	6.1446.160	<b>Thermometer stopper, B-14/15, oblique insertion</b> To have the inserts of the 6.1458.010 titration head in a sloping position.	

---

**3 pcs**

**6.2709.070**

**Guiding sleeve**

For fixation of buret and dosing tips in titration heads and electrode holders

