





Accessories for 6.2041.350

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Optional accessories

Order no.	Description	
6.1432.210	Sample beaker, clear glass, 75 mL, 1 piece	
	Sample beaker, clear glass, 75 mL, suitable for sample racks: 6.2041.340, 6.2041.350, 6.2041.800 and 6.2041.900, 1 piece	
6.1458.070	Titration head / 2 x SGJ 14	
	Titration head for 75 mL beakers. With 2 openings SGJ 14, 2 openings 6.4 mm, and 3 pre-mounted dosing tips. Suitable for rinsing and aspiration equipment.	
6.1459.400	Sample beaker, plastic (PE), 75 mL, 144 pieces	
	Sample beaker, plastic (polyethylene, PE) 75 mL, suitable for sample racks: 6.2041.340, 6.2041.350, and 6.2041.800, suitable for 6.02041.040 OMNIS Software sample rack, 144 pieces.	
6.1462.050	Robotic arm with titration head for 786 Swing Head, left or right swinging	
	Titration head for Swing Head	