





## Accessories for 2.942.0040

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery 2.942.0040

Qty.	Order no.	Description	
1 pcs	1.942.0040	<b>942 Extension Module Vario HPG</b>	 <p>The 942 Extension Module Vario HPG is an extension module for the Professional IC Vario line of instruments. It permits the addition of a further eluent in a Professional IC Vario <b>high-pressure gradient system</b>.</p> <p>Typical applications</p> <ul style="list-style-type: none"> <li>• High-pressure gradient systems with up to four eluents</li> </ul>
1 pcs	6.1602.160	<b>Bottle attachment GL 45 for eluents and reagents</b>	 <p>For eluent, reagent, and rinsing bottles, with connectors for the adsorber tube and the aspiration tubing.</p>
1 pcs	6.1608.070	<b>Eluent bottle / 2 L / GL 45</b>	
1 pcs	6.1609.000	<b>Adsorbing tube, large and bent</b>	 <p>For filling with adsorber material.</p>

---

**1 pcs**     **6.1805.060**     **FEP tubing / M6 / 60 cm**  
With light and kink protection



---

**1 pcs**     **6.1831.100**     **PEEK capillary 0.25 mm i.d., 1 m**



---

**1 pcs**     **6.1834.080**     **Aspiration tubing, 2 m**  
Aspiration tubing for Professional IC instruments



---

**1 pcs**     **6.2023.020**     **Clip for SGJ 14/15**  
Clip for SGJ 14/16



---

**1 pcs**     **6.2156.060**     **Cable Extension Module - Professional IC, 40 cm**  
Cable connecting the extension module with a Professional IC instrument



**1 pcs**     **6.2251.000**     **Colored sleeves for capillaries**  
 Colored pieces of heat shrink tubing for capillary coding. Three pieces each of five different colors.



**1 pcs**     **6.2621.030**     **Hexagon key 4 mm**



**1 pcs**     **6.2621.100**     **3 mm hex key for IC Sample Processors**



**1 pcs**     **6.2744.070**     **Pressure screw short**  
 Short version. With UNF 10/32 connection. 5 pieces. For the connection of PEEK capillaries



**1 pcs**     **6.2744.210**     **Tubing adapter for aspiration filter (ProfIC)**  
 For Professional IC instruments



**1 pcs**     **6.2758.000**     **Mixing coil for high pressure gradients**  
 For high pressure gradient mixing



---

**1 pcs**     **6.2816.020**     **Syringe 10 mL with Luer connection**  
For various applications in IC and VA



---

**1 pcs**     **6.2816.040**     **Purge tubing**  
With PTFE tubing and Luer connection. For syringes. For aspirating eluents.



---

**1 pcs**     **6.2821.090**     **Aspiration filter**  
Pore size 20 µm, set of 5 pieces. For 6.1834.000 Aspiration tubing and 6.1821.040 and 6.1821.050 Filter tubes.



---

**1 pcs**     **6.2821.130**     **Spare filter for inline filter**  
Spare filter for inline filter.



## Optional accessories

Order no.	Description	
6.2103.170	<b>Adapter to Leak Sensor at Professional IC</b> Adapter to connect two leak sensors to a Professional IC instrument	
6.2156.070	<b>Cable Extension Module - Professional IC, 1 m</b> Cable connecting the extension module with the Professional IC instrument	
6.2617.040	<b>Tool for piston seal, macro</b> For removing and assembling the piston seal for all macro pump heads	
6.2741.000	<b>Piston seal</b> For all IC pumps. Particularly well-suited for eluents with organic solvents.	
6.2741.040	<b>PE/PTFE piston seal, macro</b> For all macro pump heads	
6.2744.240	<b>Cross connector UNF 10/32 PEEK</b> Cross connector with four connectors for mixing three solutions.	

---

**6.2824.130      Macro pump head PEEK**

Macro pump head for intelligente IC instruments, flow range 0.1...20 mL/min, maximum pressure 12.5 MPa.



---

**6.5332.000      IC tubing kit for second extension module**

