










Accessories for 2.942.0500

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.942.0500

Qty.	Order no.	Description	
1 pcs	1.942.0500	942 Extension Module Vario SeS/PP	
		<p>The 942 Extension Module Vario SeS/PP is an extension module for the Professional IC Vario line of instruments. It permits the integration of a complete sequential suppression with a peristaltic pump for suppressor regeneration in a Professional IC Vario system.</p> <p>Typical applications:</p> <ul style="list-style-type: none"> • Sequential suppression for an additional analysis channel 	
1 pcs	6.1602.150	Bottle attachment / GL 45 - 3 x UNF 10/32	
		For the connection of 1/16 in. capillaries. Used in conjunction with MSM auxiliary solutions and in Metrohm Inline Sample Preparation.	
1 pcs	6.1608.020	Glass bottle / 1000 mL / GL 45	
		Bottle for auxiliary solutions.	
1 pcs	6.1803.030	PTFE capillary 0.5 mm i.d. / 3 m	
		Capillary for inline dialysis, for Dialysis Unit, IC Dialysis Sample Processor, IC Liquid Handling Dialysis Unit	

2 pcs	6.1826.420	PharMed® pump tubing (orange/yellow), 3-stopper	
For MSM regeneration with peristaltic pump.			
1 pcs	6.1831.070	PEEK capillary 0.75 mm i.d., 70 cm	
For Advanced VA Sample Processor			
1 pcs	6.2156.070	Cable Extension Module - Professional IC, 1 m	
Cable connecting the extension module with the Professional IC instrument			
1 pcs	6.2251.000	Colored sleeves for capillaries	
Colored pieces of heat shrink tubing for capillary coding. Three pieces each of five different colors.			
1 pcs	6.2744.034	Coupling olive/UNF 10/32, 2x	
Connection between pressure screw and pump tubing. 2 pieces. For IC instruments with peristaltic pump.			
2 pcs	6.2744.180	Pump tubing connector with security lock and filter	
For the connection of a pump tubing and a capillary with built-in filter			
1 pcs	6.5000.020	Accessory kit: Vario/Flex SeS	
Set of accessories for sequential suppression (SeS) for Professional IC Vario and Compact IC Flex instruments. Without Suppressor Rotor.			

Optional accessories

Order no.	Description	
6.2103.170	Adapter to Leak Sensor at Professional IC	
	Adapter to connect two leak sensors to a Professional IC instrument	
6.2832.000	MSM Rotor A	
	Suppressor rotor for all IC instruments with MSM (Metrohm Suppressor Module)	
6.2842.000	MSM-HC Rotor A	
	Suppressor rotor for all IC instruments with MSM-HC (Metrohm Suppressor Module with high capacity)	
6.2842.020	Adapter sleeve for Suppressor Vario	
	The use of the adapter sleeve makes it possible to install both the Standard MSM and the MSM-LC rotor in the MSM-HC housing.	
6.2842.200	MSM-HC Rotor C	
	Suppressor rotor for all IC instruments with MSM-HC drive (Metrohm Suppressor Module with high capacity) for cation suppression.	

6.2844.000

MSM-LC Rotor A

Suppressor rotor for use with 2 mm columns.

