


## Accessories for SPELECNIR

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery SPELECNIR

Qty.	Order no.	Description	
1 pcs	DropView SPELEC	DropView SPELEC Software	
		DropView SPELEC is a Spectroelectrochemical software that controls SPELEC instrument, offering a perfect synchronization of the optical and electrochemical measurements, as well as advanced tools for data treatment.	

---