




## Accessories for 2.700.0010

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery 2.700.0010

Qty.	Order no.	Description	
1 pcs	1.700.0010	<b>700 Dosino</b> Drive for Dosing Unit, with fixed cable and 9-pin plug (DB9).	
1 pcs	6.1546.030	<b>Extractor rod</b> For PTFE piston of dosing unit	
1 pcs	6.2134.200	<b>Adapter DB9 St. /Mini Din8 Bu</b>	
1 pcs	6.2148.040	<b>Adapter for Dosino 700 (series 23 or higher)</b>	

## Optional accessories

Order no.	Description
<b>6.2134.020</b>	<b>Adapter DB9 female/mini-DIN8 male</b>
DB9 socket and mini-DIN 8 plug. Connection 700.0010 Dosino to 736, 751, 758, 799 Titrinos, 756, 831 KF Coulometer, 711 Liquino.	

