











Accessories for 2.765.0040

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.765.0040

Qty.	Order no.	Description	
1 pcs	1.728.0010	728 Stirrer without stand For mounting on a stand rod 10 mm dia.; without base plate (or holder) and without stand rod; without power adapter, with fixed connecting cable.	
1 pcs	1.765.0020	765 Dosimat Universal dispensing unit for liquid handling in laboratories. Including push-button cable for manual dispensing control, keypad and 6.3026.220 Exchange Unit 20 mL. Dosimat with analog output for volume registration by recorders.	
1 pcs	6.1903.020	Stirring bar / 16 mm Stirring Bar with magnetic core, PTFE coated, length 16 mm	
1 pcs	6.1903.030	Stirring bar / 25 mm Stirring bar with magnetic core, PTFE covering, length 25 mm.	

1 pcs	6.2001.010	Support stand with support rod For Dosimats and Titrinos	
1 pcs	6.2013.010	Clamping ring For support rods with a diameter of 10 mm.	
1 pcs	6.2021.020	Electrode holder Electrode holder for 4 electrodes and 2 buret tips	
1 pcs	6.2107.000	Push-button cable for Dosimat For manual operation of Dosimats	
1 pcs	6.2149.000	Keyboard for 765/776 Keyboard for 765, 776 Dosimat	
1 pcs	6.2723.130	Cover for 655/665/765 Dust cover for 655, 665, 765 Dosimats	

1 pcs

6.2739.010

Wrench

For exchange units.



1 pcs

6.3026.220

Exchange Unit 20 mL

Exchange unit with integrated data chip with 20 mL glass cylinder and light protection. PCTFE/PTFE flat stopcock, FEP tubing connection, anti-diffusion buret tip and standard amber glass reagent bottle.

