





## Accessories for 2.768.0010

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery 2.768.0010

Qty.	Order no.	Description	
1 pcs	6.1448.040 Set of 5 pcs.	Septum 22 mm, 5 pieces	
1 pcs	6.1543.060	Tip / M6 thread  Tip with M6 thread. Together with the 6.1805.160 capillary and the 6.1446.030 link stopper, it makes up the entire 6.1537.010 buret tip.	
2 pcs	6.1602.145	Drying flask cap  Screw cap with GL 45 thread for the 6.1608.050 drying flask.	
2 pcs	6.1608.050	Drying bottle / 100 mL / GL 45	

---

**1 pcs**      **6.1805.040**      **FEP tubing / M6 / 9 cm**  
With light and kink protection



---

**1 pcs**      **6.1805.070**      **FEP tubing / M6 / 10.5 cm**  
With light and kink protection



---

**2 pcs**      **6.1805.080**      **FEP tubing / M6 / 25 cm**  
With light and kink protection.



---

**3 pcs**      **6.1805.180**      **FEP tubing / M6 / 16 cm**  
With light and kink protection



---

**1 pcs**     **6.1808.020**     **Tubing adapter olive / M6 inner**  
1 M6 inner thread and 1 tubing nozzle. Part of the screw connections for Exchange Units and stability measuring instruments.



---

**2 pcs**     **6.1821.040**     **Filter tube**  
Filter tube for 6.1608.050 Drying bottle. For Rancimats and Karl Fischer ovens.



---

**1 pcs**     **6.2104.080**     **Electrode cable 1 m, 2xB**  
With banana plug 4 mm dia. For connecting electrodes with Metrohm plug-in head G to measuring instrument.



---

**1 pcs**     **6.2701.060**     **GL 25 screw cap with hole**  
Made of PBT. With a hole for the 6.1448.040 septum. For 707, 768 KF Ovens



---

**2 pcs**     **6.2724.010**     **Dust filter**  
Dust filter for Rancimats and Karl Fischer ovens.



---

**1 pcs**

**6.2730.050**

**Stopper**

With E.301.0041 O-ring. For using 6.1824.010 Double buret tip in 6.1414.030 KF titration vessel lid



---

**1 pcs**






**6.2811.000**

**Molecular sieve**

Molecular sieve. Bottle containing 250 g. Pore size: 0.3 nm. Without moisture indicator. For Rancimats and Karl Fischer instruments.



## Optional accessories

Order no.	Description	
<b>6.1446.170</b>	<b>Heating tube stopper</b> for coulometric KF vessels in connection with heatable tubing connection	
<b>6.2125.010</b>	<b>RS-232 cable 3 m, 25-pin – 9-pin</b> Connecting RS-232 interface to IBM PC AT or compatible computer. Adapter 25-pin to 9-pin	
<b>6.2125.050</b>	<b>RS-232 cable, 25-pin plug</b> Connecting Metrohm RS-232 interface to Epson printer LX 300. Used with 691, 692, 713 pH Meters, Titrinos, and 737 KF Coulometer.	
<b>6.2125.060</b>	<b>Connection cable RS-232 to IBM PC</b> Connecting Metrohm RS-232 interface to PC. Used with 691, 692, 713 pH Meters, Titrinos, and 737 KF Coulometer	
<b>6.2125.100</b>	<b>RS-232 cable 2 m, 25-pin D – sub 25</b> RS-232 extension cable (1:1) with special EMC shielding.	
<b>6.2125.110</b>	<b>RS-232 cable 3 m, 25-pin – 9-pin</b> Connecting Metrohm RS-232 interface to PC. Used with pH Meters, Titrinos, KF Coulometers, VA Trace Analyzer.	
<b>6.2125.130</b>	<b>RS-232 cable, 1.8 m, DB25 (socket) – DB9 (plug)</b> Connecting Metrohm RS-232 interface (DB25) to a Custom DP40 printer or Mettler PG 5200 balance.	

---

**6.2141.010**      **Cable 787/795/758/831/756 - 707/768**

Connecting 707, 768 KF Oven to KF Titrinos or 756 KF Coulometer. For older KF Coulometers see 6.2141.000

