

Accessories for 2.772.0120

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.772.0120

| Qty. | Order no. | Description | |
|-------|------------|---|---|
| 1 pcs | 1.772.0110 | 772 Pump Unit Peristaltic pump for direct connection to Sample Processors. |  |
| 1 pcs | 6.1446.160 | Thermometer stopper, B-14/15, oblique insertion To have the inserts of the 6.1458.010 titration head in a sloping position. |  |
| 1 pcs | 6.1543.170 | Aspiration tip / M8 thread Aspiration tip for the Sample Processors. |  |
| 1 pcs | 6.1621.000 | Container 10 L As rinsing or waste container in automated systems. |  |

2 pcs **6.1801.120** **PVC tubing / 4 mm / 6 mm / 2 m**



1 pcs **6.1805.510** **PTFE tubing / M8 / 60 cm**

With protection against kinks.



1 pcs **6.1812.000** **PTFE tubing 4 mm / 6 mm / 4 m**



2 pcs **6.1820.050** **M6 / Olive adapter**

For 6.1826.100 pump tubing, to peristaltic pumps.



1 pcs **6.1820.070** **Coupling 4 mm / 6 mm / M8 inner**

For extending aspiration tubing in automated systems. Connects the M8 aspiration tubing directly to the 6.1812.000 PTFE tubing.



1 pcs **6.1826.100** **PP pump tubing 6.4 mm / 1 m**

Replaces the pump tubing used with the 772 Pump Unit and the 843 Pump Station (peristaltic).





1 pcs **6.1828.000** **PVDF connection nipple**

For 6.1621.000 canister.



Optional accessories

| Order no. | Description | |
|--|---|---|
| 6.1826.160 | Viton pump tubing |  |
| For peristaltic pumps. | | |
| 6.2160.010 | Adapter cable / 2 m / M8 socket / 2 x plug B |  |
| Adapter cable / 2 m / M8 socket / 2 x plug B. For 772 Pump Unit. | | |