

Accessories for 2.812.0020

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.812.0020

Qty.	Order no.	Description	
1 pcs	1.812.0020	812 IC Valve Unit – 2 metal-free 6-way injection valves Instrument with two metal-free, electrically operated six-way injection valves. Particularly suitable for sample injection, column switching and extension to an inline dialysis system, for example. Simple remote control via Remote lines.	
1 pcs	6.1803.020	PTFE capillary 0.97 mm i.d. / 5 m For all IC instruments	
1 pcs	6.1825.210	PEEK sample loop 20 µL For injection valve, with 2 PEEK pressure screws	
1 pcs	6.1831.010	PEEK capillary 0.25 mm i.d. / 3 m For all IC components.	

2 pcs **6.2744.010** **Pressure screw 5x**
With UNF 10/32 connection. For the connection of PEEK capillaries



1 pcs **6.2744.120** **Luer F pressure screw adapter**
Connection between pressure screw and syringe.



1 pcs **6.2816.020** **Syringe 10 mL with Luer connection**
For various applications in IC and VA



Optional accessories

Order no.	Description	
6.1803.030	PTFE capillary 0.5 mm i.d. / 3 m Capillary for inline dialysis, for Dialysis Unit, IC Dialysis Sample Processor, IC Liquid Handling Dialysis Unit	
6.1825.220	PEEK sample loop 100 µL For injection valve, with 2 PEEK pressure screws	
6.1825.230	PEEK sample loop 10 µL For injection valve, with 2 PEEK pressure screws	
6.5904.030	Rotor for flat plate injector (metal-free) For all IC instruments with metal-free injectors	
6.5904.050	Injector with 0.25 µL internal loop Flat plate injector with 0.25 µL internal loop	

6.5904.060 **Rotor (internal loop, 0.25 µL)**

Rotor for 4-port flat plate injector (6.5904.050), 0.25 µL internal loop



6.5904.070 **Rotor (internal loop, 1.0 µL)**

Rotor for 4-port flat plate injector (6.9956.001), 1.0 µL internal loop

